



Breathable
Membrane



Easy to Use



Leak Proof




Antibacterial



Durable
Long Life



PVC Medical
Grade Material



AM-101 Air Mattress

Features:

- Reliable & Easy to use.
- Low noise and low vibrations.
- Knob to adjust pressure.
- Hooks to hand the motor easily on the bed frame.

Specifications:

WEIGHT CAPACITY:

Maximum Patient Weight

130kg.

ELECTRICAL:

Operating Power

~220V. 50Hz. 0.1A (Max)

ENVIRONMENTAL CONDITIONS:

Ambient Temperature

10 °C to 40 °C

Relative Humidity

10% to 75%

STORAGE AND SHIPPING CONDITIONS:

Ambient Temperature

~18 °C to 43 °C

Relative Humidity

10% to 795%

SPECIFICATION:

Cycle Time

6 Minutes

Air Output

6-7 L/Minutes

Pressure Range

70-130 mmHg

How to Use:



1. Unfold & lay the bubble mattress on a foam mattress in a way that the air connections are on the foot side.



2. Place the pump on a flat table or hang it on the bed with the help of the hanger at the foot side.



3. Connect one end of the both air tubes with the air pump and the other end with the mattress. Connect them firmly at the ends. Make sure that air pipes are not cinked.



4. Switch on the electric pump by connecting with 220-V AC. The mattress will start inflating. Keep the pressure adjusting knob at full in the beginning. This will take around 10 minutes to inflate.



5. Adjust pressure ideally with the adjusting Knob in such a way that there is a two fingers gap between the deflated air bubbles and the patient.

Contact us:



Mobile No./ WhatsApp No. : +91 9310113243



Web : www.klifecare.com



Email : info@klifecare.com

Kannu Impex (India) Pvt. Ltd.



412, Suncity Trade Tower, Sector-21, Gurugram-122015