



Best Quality



Leak Proof



Fast Relief



Easy to Use

HWB-101 Hot Water Bag

Features:

- Treats Muscular Pain
- Induces Relaxation
- Releases Stress Instantly
- Relieves Joint Pains
- Muscular Cramps
- Menstrual Pains
- Stomach and Back Aches

Specifications:


Material	Natural Rubber
Temperature	100°C to 42°C
Using Time	15 to 30 Minutes

Warning

Hot water bag can be dangerous & cause burns if they are not cared for or used in the proper way. Please read our guides for filling, using and caring for your hot water bag below.

How to Use:

1



When filling your hot water bag, do not use boiling water. Once the water has boiled allow it to stand for a couple of minutes before filling your bag. Filling a bag with boiling water can cause splash back, which may cause burns.

2



When filling your bag hold the bag by the neck in an upright position and fill slowly to avoid water splashing back.

3




Only fill the bag to the MAXIMUM of two-thirds capacity.

4



ALWAYS expel air from the bag by lowering it carefully onto a flat surface until the water appears at the opening.

5



Screw the stopper sufficiently tight to ensure there is no leakage. Finger tight should be adequate.

6

Hot water bag should always be filled by a responsible adult.

Contact us:



Mobile No./ WhatsApp No. : +91 9310113243




Web : www.klifecare.com



Email : info@klifecare.com

Kannu Impex (India) Pvt. Ltd.

 412, Suncity Trade Tower, Sector-21, Gurugram-122015