

Doc No – Ab/TDS/CT-01

Rev No: 08 Date- 02/2022-23

TECHNICAL DATA SHEET

Abdos Labtech Pvt. Ltd.

Micro Centrifuge Tube, PP

Valid for Cat. No. P10201, P10202, P10203, P10204, P10204A, P10207, P10201A, P10205, P10203A

1. Description/Specification

1.1 Description:

Micro Centrifuge tube PP, DNase, RNase, Human DNA, ATP & Pyrogen Free , CE & IVD certified

1.2 Design:

- 1.2.1 Tube: Cylindrical upper and conical bottom
- 1.2.2 Tube cap: Attached Snap Cap with lid open force kgf

1.3 Performance Specification:

Code	Capacity (ml)	OD (mm)	Height (mm)	Sterile	Colour	No./Case
P10201	0.5	10	31.5	NO	Natural	1000
P10202	1.5	13	41	NO	Natural	500
P10203	2.0	13	41	NO	Natural	500
P10207	0.5	10	31.5	NO	Red, Blue, Yellow & Green	1000
P10201A	0.5	10	31.5	NO	Amber	1000
P10204	1.5	13	41	NO	Red, Blue, Yellow & Green	500
P10204A	1.5	13	41	NO	Amber	500
P10205	2.0	13	41	NO	Red, Blue, Yellow & Green	500
P10205A	2.0	13	41	NO	Amber	500

1.4 Material:

- 1.4.1 Polypropylene (PP) confirm US FDA 21CFR
- 1.4.2 Complies with the requirement of European Pharmacopoeia, Plastic Containers and closures for Pharmaceutical use.

Ph. No: +91- (0) 1332235277 Email:labtech@abdosindia.com Web: www.abdoslifesiences.com GST: 05AAGCA8605A1Z7

www.rajbiosis.com +91-998-3333-469



- 1.4.3 The product has passed the USP testing including class VI tests. Materials has successful passed the biological tests according to ISO 10993 external communication devices for indirect blood contact for Prolonged period.
- 1.4.4 Free from Heavy metals, RoHS Free according to meets relevant requirements of EU regulations1935/2004/CE ,2002/95/EC as amended, UE 10/2011 as amended. EU CONEG & it's following amendments. It is safe use in contact with food stuff, Pharmaceuticals & drinking water
- 1.4.5 Product is manufactured without the use of Raw Materials of animal origin. Therefore, BSE/TSE issue does not concern this product to the best of our Knowledge. No materials containing latex or natural rubber are used in the Manufacturing, handling & packaging process for this product.

1.5 Color:

Translucent (Natural), Coloured (Red, Blue, Yellow & Green) and Amber

1.6 Graduation and Writing Surface

- 1.6.1 Large Molded Graduation
- 1.6.2 Frosted Surface for easy writing on tube

1.7 Quality control:

Product are manufactured, inspected and found to be in compliance with product & quality specification requirements as documented in ISO 9001:2015 & 13485:2016 Quality Management System. Testing of various characteristics in accordance with the valid specification, inspection reviewed & signed for product release.

1.8 Other information

Recommended for Single use only

1.9 Basic Features

1.9.1 Temperature range: Tube:-80 0 C to + 121 0 C

1.9.2 Autoclavability: Yes, Autoclavable at 15 psi and 121 °C for 20 minutes

1.9.3 Centrifugation: 14,000* to 20,000 RCF (* Rated with Phenol/ Chloroform)

1.9.4 Leak Test at 635 mm Hg: No leakage

1.9.5 Boiling test: Pass at 100°C for 30 minutes in water bath.

1.9.6 Manufactured in Clean room ISO Class 8 (100,000)

1.9.8 Chemical resistance: see website (www.abdoslifesciences.com)

Autoclaving Guideline for Plastics

All plastic consumable is Autoclavable at 121°C for 20 minutes, 15 psi/1atm. Do not use dry cycle as this may cause deformation or weakening of plastic. Allow temperature within the autoclave to return to atleast 50°C before removing product.

Ph. No: +91- (0) 1332235277 Email:labtech@abdosindia.com Web: www.abdoslifesiences.com GST: 05AAGCA8605A1Z7



Instruction for Use

Product can be used for food contact also.

Storage temp. Max 30°C

Before using the product check crackness, pouch punctured etc.

Check all the information printed on pouches or boxes.

Always wear gloves while handling specimen.

Ensure all sample tubes are evenly filled


For each tube inserted in the rotor, add a tube of equal weight directly opposite it.

WARNING

DO NOT store strong oxidizing agents.

DO NOT place any plastic labware in a flame.

DO NOT mix any chemicals that may result in a thermic reaction, which can cause product failure Wear a face shield and / or safety goggles if you have to work anywhere near a centrifuge that's in use. **Do not bump, jar, or move the centrifuge while the rotor is spinning**. Make sure you don't have the cord dangling from a table edge where someone could catch their foot in it and pull down the centrifuge Consult local fire codes prior to storage of flammable liquids.



www.rajbiosis.com +91-998-3333-469



Abdos Labtech Private Limited
(A Divison of Abdos Group)
(An ISO 9001,13485,14001,18385,GMP & CE certified company)
Shiv Ganga Industrial Estate
Plot No-1, Lakheswari, Bhagwanpur,
Roorkee-247667, (UK) INDIA

Registered office
13A Government
Place East, 2nd Floor Kolkata
West Bengal – 700069
CIN: U85100WB2008PTC123903
Corporate Office
319, 3rd Floor DIE Tower A DDA Dies

319, 3rd Floor, DLF Tower A, DDA Dist. Centre Jasola, New Delhi - 110025

Ph. No: +91- (0) 1332235277 Email:labtech@abdosindia.com Web: www.abdoslifesiences.com GST: 05AAGCA8605A1Z7





www.rajbiosis.com +91-998-3333-469



319, 3rd Floor, DLF Tower A, DDA Dist. Centre Jasola, New Delhi - 110025 Email:labtech@abdosindia.com Web: www.abdoslifesiences.com GST: 05AAGCA8605A1Z7

Ph. No: +91- (0) 1332235277