

DISINTEGRATION TEST APPARTUS MODEL : VTS-P



This disintegration tester determines whether suppositories soften or disintegrate within a prescribed time when placed in an immersion fluid under the experimental conditions.

- Designed for the measurement of disintegration time of pessaries and suppositories.
- LCD display screen to indicate test parameters.
- Automatic rotation through 180 deg and selectable cycle time.
- Built-in circulation pump and sensitive circuit to maintain accurate temperature between 36.0° - 37.0°C.
- Three samples are tested at the same time.
- The instrument meets requirements of pharmacopoeias.

Attachment to test Disintegration Test time of pessaries.



Stand For Suppository



Stand For Suppository



Baskets For Pessaries



Baskets For Suppository
& Pessaries

VEEGO INSTRUMENTS CORPORATION

201/202, Unique Industrial Estate, Chakala Road, Andheri (E), Mumbai - 400099, India
Call : +91 22 4023 2455 / 4010 3727 | Email : office@veego.in | www.veego.in