

A complete CSSD solution for your hospital



AUTOCLAVE

Vetrical Autoclave Cylindrical Horizontal Autoclave Rectangular Horizontal Autoclave Sliding Horizontal Autoclave Bio MEDICAL Waste Autoclave Industrial Autoclave Table TOP Autoclave

ETO STERILIZER

Ethylene Oxide Gas Sterilizer (ETO STERILIZER)

Industrial (ETO STERILIZER) Table Top (ETO STERILIZER)

PRODUCTS

Ultra Sonic Clener Mortuary Chamber LaminorBOD

LAUNDRY STATION

Washing Machine Hydro Extractor Tumble Dryer Hydro Machine Bhagona

COMPLETE CSSD SOLUTION FOR YOUR HOSPITAL



WELCOME TO YHE WOELD OF - CSSD -

CATALOG CONTAINS

- All Vertical Autoclave
- Horizontal Autoclave
- ETO Sterilizer
- Semi Automatic
- Fully Automatic
- Single Door, Double Door

COMPLETE SETUP OF LOUNDRY STATION ON ROOF

- Washing Machine
- Hydro Extractor
- Tumble Dryer
- Hydro Machine
- Vacuum Press
- etc...

Vertical Triple Walled High Pressure Autoclave Model Number: NEXX-001

NEXXGEN SURGICAL, Have Tow Model of Vertical Triple Walled Autoclave Will Be TOP Loading and in Vertical Shape With Triple Layer With Boiler. These High Pressure TOP Loading Autoclaves Play Important Role in Sterilization Procedures Required in Hospitals, Pharmaceutical Companies, Biopharma and Microbiology Laboratories etc.

NEXXGEN SURGICAL Vertical Autoclave Are Made of Stainless Steel to Ensure Corrosion Resistance Surface. Tipple Walled Construction Ensures Excellent Level of Safety. Boiler is Fitted With Chamber and Connected to the Chamber. Water Level Indicator is Also Attached to the Boiler With Automatic Low Water Cut-off.

User Can Easily Set The Desired Temperature With The Help Of Digital Controller Fitted on The Front Side of The Unit. We Provide Radial Locking Pressure Door With Silicon Gas Kit. We Use Automatic Pressure Control Switch an User Easly Can Set Temperature by Digital Controller. Our Machine Work With Pressure Between 15-22 PSI.

FEATURES & BENEFITS:

- Vertical Triple Walled Autoclave
- Stainless Steel Interior and Exterior
- Triple Walled Construction
- Standard Temperature 121°C
- Working Pressure 15 psi to 22 psi
- Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch

TECHNICAL SPECIFICATIONS

Model	NEXX-001
Design	Vertical
Construction	Triple Walled
Working Pressure	Variable 15 to 22 PSI
Working Temperature	121ºC
Internal Chamber	Corrosion Resistant Stainless Steel SS 316
External Chamber	Corrosion Resistant 304 Grade Stainless Steel
Finishing	Mirror Finished Inner and External Surface
All Piping	SS 304
Door	Single Doors
Sterilization Period	25 to 35 Minutes
Steam Exhaust	5 to 7 Minutes
Boiler	 Made of SS 304 Sheet Water Inlet & Outlet Valves Water Level Glass Window Boiler Safety Valve With ISI Marked Water Immersion Type Industrial Heating Elements.
Hydraulic Tested	on 40 psi
Sensor	PT 100
Safety Valve	For Excess Pressure Inside Chamber and Steam Generator.
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes
Optional Accessories	 PLC Operation Vacuum Drying Printer Auto Add Water Door Inter Lock Automatic Pressure Control Switch Water Cut - off Temperature Indicator Caster Wheels
Power Supply	220 Volts / 440 Volts (Three Phase)

- Vacuum Breaker in Chamber
- ISI Safety Valves
- Digital Temperature Controller with Timer and Alarm (Buzzer)
- Automatic Water Level Cut-off
- Water Level Glass To Check Water Level
- SS Pressure Gauges
- ISI Marked Flange / Immersion Heating Elements



NEXX-001

Size of Vertical Autoclave

ne and terilization	INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
	250 X 450 MM	10" X 18"	22	2
	300 X 500 MM	12" X 20"	40	3
	350 X 550 MM	14" X 22"	50	3
ntrol Switch	400 X 600 MM	16" X 24"	78	4
	450 X 600 MM	18" X 24"	98	4
ee Phase)	550 X 750 MM	22" X 30"	180	6

Steam Storage High Pressure Vertical Autoclave Model Number: NEXX-001A

NEXXGEN SURGICAL, Have Tow Model of Vertical Triple Walled Steam Storage Autoclave Will Be TOP Loading and in Vertical Shape With Two Chamber. One is Inner Chamber For Equipmets Sterilization and Second Will Store Steam in Jacket With Boiler. These High Pressure TOP Loading Autoclaves Play Important Role in Sterilization Procedures Required in Hospitals, Pharmaceutical Companies, Biopharma and Microbiology Laboratories etc.

NEXXGEN SURGICAL Vertical Autoclave Are Made of Stainless Steel to Ensure Corrosion Resistance Surface. Tipple Walled Construction Ensures Excellent Level of Safety. Boiler is Fitted With Chamber and Connected to the Chamber. Water Level Indicator is Also Attached to the Boiler With Automatic Low Water Cut-off.

User Can Easily Set The Desired Temperature With The Help Of Digital Controller Fitted on The Front Side of The Unit. We Provide Radial Locking Pressure Door With Silicon Gas Kit. We Use Automatic Pressure Control Switch an User Easly Can Set Temperature by Digital Controller. Our Machine Work With Pressure Between 15-22 PSI.

FEATURES & BENEFITS:

- Vertical Triple Walled Steam Storage Autoclave
- Stainless Steel Interior and Exterior
- Triple Walled Construction
- Standard Temperature 121°C
- Working Pressure 15 psi to 22 psi
- Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch

TECHNICAL SPECIFICATIONS

Model	NEXX-001A
Design	Vertical
Construction	Triple Walled
Working Pressure	Variable 15 to 22 PSI
Working Temperature	121ºC
Internal Chamber	Corrosion Resistant Stainless Steel SS 316
External Chamber	Corrosion Resistant 304 Grade Stainless Steel
Finishing	Mirror Finished Inner and External Surface
All Piping	SS 304
Door	Single Doors
Sterilization Period	25 to 35 Minutes
Steam Exhaust	5 to 7 Minutes
Boiler	 Made of SS 304 Sheet Water Inlet & Outlet Valves Water Level Glass Window Boiler Safety Valve With ISI Marked Water Immersion Type Industrial Heating Elements.
Hydraulic Tested	on 40 psi
Sensor	PT 100
Safety Valve	For Excess Pressure Inside Chamber and Steam Generator.
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes
Optional Accessories	 PLC Operation Vacuum Drying Printer Auto Add Water Door Inter Lock Automatic Pressure Control Switch Water Cut - off Temperature Indicator Caster Wheels
Power Supply	220 Volts / 440 Volts (Three Phase)

- Vacuum Breaker in Chamber
- ISI Safety Valves
- Digital Temperature Controller with Timer and Alarm (Buzzer)
- Automatic Water Level Cut-off
- Water Level Glass To Check Water Level
- SS Pressure Gauges
- ISI Marked Flange / Immersion Heating Elements



NEXX-001A

tion	INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
	250 X 450 MM	10" X 18"	22	2
	300 X 500 MM	12" X 20"	40	3
	350 X 550 MM	14" X 22"	50	3
witch	400 X 600 MM	16" X 24"	78	4
	450 X 600 MM	18" X 24"	98	4
se)	550 X 750 MM	22" X 30"	180	6

Cylindrical Horizontal Semi Automatic Autoclave Model Number: NEXX-002

NEXXGEN SURGICAL, Have Comprehensive Range Of Horizontal Autoclaves. Cylindrical Horizontal Autoclave Will Be Front Loading And In Horizontal Shape With Triple Layer With Boiler. These High Pressure Front Loading Autoclaves Play Important Role In Sterilization Procedures Required In Hospitals, Pharmaceutical Companies, Biopharma And Microbiology Laboratories Etc.

NEXXGEN SURGICAL Horizontal Autoclave Units Are Made Of Stainless Steel To Ensure Corrosion Resistance Surface. Tipple Walled Construction Ensures Excellent Level Of Safety. Boiler Is Fitted Separately And Connected To The Chamber. Water Level Indicator Is Also Attached To The Boiler With Automatic Low Water Cut-off.

User can easily set the desired temperature with the help of digital controller fitted on the front side of the unit. we provide radial locking pressure door with silicon gas kit. we use automatic pressure control switch and user easly can set temperature by Digital Controller. our machine work with pressure between 15-34 PSI.

TECHNICAL SPECIFICATIONS

Model	NEXX-002	
Design	Cylindrical Horizontal Type	
Construction	Triple Walled	
Valve for Jacket to Chamber Transfer, Chamber Release and For Post Vacuum or Dary Mode	Use Heavy Duty Multi-port Valve	
Self Steam Generator System	SS 304 Boiler With Heater	
Pre-Vacuum or High Vacuum	Yes	
Water in Boiler	Distilled / Normal Water	
Heater	Industrial Immersion Type Water Heater	
Working Pressure	Variable 15 to 34 PSI	
Working Temperature	121°C (up Graded up to 135°C	
Internal Chamber	Corrosion Resistant Stainless Steel SS 316	
External Chamber	Corrosion Resistant 304 Grade Stainless Steel	
Finishing	Mirror Finished Inner and External Surface SS 304 / 316L	
All Piping	SS 316L	
Door	Single or Double Doors SS 304	
Insulation	Glass Wool Fiber 100mm	
Sterilization Period	25 to 30 Minutes	
Steam Exhaust	5 to 7 Minutes	
	Decontaminated by Means of a HEPA Filter	
Air Removed During The Vacuum Cycle	(Class H13 or Higher)	
Dry	5 - 10 Minutes	
Boiler	 Made of SS 304 Sheet Water Inlet & Outlet Valves Water Level Glass Window Boiler Safty Valve With ISI Marked Water Immersion Type Industrial Heating Elements. 	
Chamber Condensate Line	Fitted With Steam Trap & NRV Value	
Hydraulic Tested	On 50 psi	
Sensor	PT 100	
Temperature	Digital Dispaly	
Temperature indication	+1 % Accuracy	
Jacket Pressure	Use Pressure Gauge	
Chamber Pressure (for Positive Pressure)	Use Compound Gauge (for Positive & Negative Pressure)	
Pressure Indication	+1.6 % Accuracy	
Safety Value	For Excess Pressure Inside Chamber, Jacket and Steam Generator.	
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Steilization Cycle Completes	
	 SS 304 Sterilizing Tray x 1pcs 	
	Water Level Glass Tube x 2pcs	
Standard Accessories	Heating Element x 1 pc	
	Silicone Door Seal x 1pc	
	PLC Operation	
	Vacuum Drying Printer	
Optional Accessories	Auto Add Water Door Interlock	
	 Automatic Pressure Control Switch Water Cut-off 	
	 Tomporatura Indiantar 	
	Temperature IndicatorCaster Wheels	

FEATURES & BENEFITS:

- Cylindrical Horizontal Autoclave
- Stainless steel interior and exterior
- Triple walled construction
- Self-locking safety door cannot be opened when chamber is under pressure
- Standard temperature 121°C to 134°C
- Working pressure 15 psi to 34 psi
- Silicon jointless rubber door gasket
- Automatic pressure control switch
- Vacuum breaker Steam trap & NRV ISI safety valves
- Digital Temperature controller with Timer and alarm (buzzer)
- Automatic water level cut-off
- Water level Glass to Check Water Level
- SS Pressure & Compound gauges
- ISI marked flange / immersion heating elements



NEXX-002

ıre)	Size of Vertical Autoclave			
namber,	INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
nd Activates e Completes	400 x 600 MM	16" X 24"	78	6-8
cs	400 X 900 MM	16" X 36"	120	8-10
DCS	400 X 1200 MM	16" X 48"	165	9
	500 X 600 MM	20" X 24"	120	8
	500 X 900 MM	20" X 36"	165	9
Switch	500 X 1200 MM	20" X 48"	250	12
ase)	550 X 1300 MM		310	

Fully Automatic LCD Controller Based Fully Covered Autoclave Model No.: NEXX-LCD01

NEXXGEN SURGICAL, Have Comprehensive Range of Horizontal Autoclaves. Fully Automatic Cylindrical Horizontal Autoclave Will Be Front Loading and in Horizobtal Shape with Layer With Boiler.

These High Pressure Front Loading Autoclaves Play Important Role in Sterilization Procedures Required in Hospitals, Pharmaceutical Companies, Biopharma and Microbiology Laboratories.

NEXXGEN SURGICAL Horizontal Autoclave Units Are Made Of Stainless Steel To Ensure Corrosion Resistance Surface Tipple Walled Construction Ensures Excellent Level Of Safety. Boiler Is Fitted Separately And Connected To The Chamber Water Level Indicator Is Also Attached To The Boiler With Automatic Low Water Cut-off

User Can Easily Set The Desired Temperature With The Help Of Digital Controller Fitted On The Front Side Of The Unit We Provide Radial Locking Pressure Door With Silicon Gas Kit We Use Automatic Pressure Control Switch And User Easly Can Set Temperature By Digital Controller Our Machine Work With Pressure Between 15-34 Psi.

FEATURES & BENEFITS:

- Fully Automatic Cylindrical Horizontal Autoclave
- FULLY COVERED AUTOCLAVE
- LCD BASED MICROPROCESSOR CONTROLLER
- Thermal Printer
- ELECTRIC Solenoid Valve
- HEAVY Vaccum EJECTOR
- Water Filling Pump
- Pre Vaccum
- Post Vaccum
- Stainless steel interior and exterior
- Triple walled construction
- Steam Stored in Jacject
- Inner Chamber made by Ss316

- Jacked & Outer Made By Ss304
- All Piping Fitting Ss304
- Standard Temperature 121°c To 134°c
- Working Pressure 15 Psi To 34 Psi
- Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch
- Vacuum Breaker
- Steam Trap & NRV
- ISI Safety Valves
- Digital Temperature Controller With Timer And Alarm (buzzer)
- Automatic Water Level Cut-off
- Water Level Glass To Check Water Level
- SS Pressure & Compound Gauges
- IS1 Marked Flange / Immersion Heating Elements





	SIZE OF INNER CHAMBER	CAPACITY	LOAD
	500 X 900MM	165LTR	12KW
001	500 X 1200MM	220LTR	12KW
	550 X 1200MM	285LTR	18KW
	600 X 1000MM	285LTR	18KW
	600 X 1100MM	320LTR	18KW

Cylindrical Horizontal Fully Automatic Autoclave Model Number: NEXX-002A

NEXXGEN SURGICAL, Have Comprehensive Range Of Horizontal Autoclaves. Cylindrical Horizontal Autoclave Will Be Front Loading And In Horizontal Shape With Triple Layer With Boiler. These High Pressure Front Loading Autoclaves Play Important Role In Sterilization Procedures Required In Hospitals, Pharmaceutical Companies, Biopharma And Microbiology Laboratories Etc.

NEXXGEN SURGICAL Horizontal Autoclave Units Are Made Of Stainless Steel To Ensure Corrosion Resistance Surface. Tipple Walled Construction Ensures Excellent Level Of Safety. Boiler Is Fitted Separately And Connected To The Chamber. Water Level Indicator Is Also Attached To The Boiler With Automatic Low Water Cut-off.

User Can Easily Set The Desired Temperature With The Help Of Digital Controller Fitted On The Front Side Of The Unit. We Provide Radial Locking Pressure Door With Silicon Gas Kit. We Use Automatic Pressure Control Switch And User Easly Can Set Temperature By Digital Controller. Our Machine Work With Pressure Between 15-34 PSI.

Our Horizontal autoclaves chamber and boiler are Hydraulic tested on 50 PSI

TECHNICAL SPECIFICATIONS

Model	NEXX-002A		
Design	Cylindrical Horizontal Type		
Construction	Triple Walled		
Automation Hardware	PLC / HMI		
Display	7"(inches) Touch Screen Multi-colour. The Information Shows In Display Like Cycle Stage, Chamber Temperature, Chamber Pressure, Jacket Pressure, Along With The Information About Failures And Interrupts.		
Printer	Thermal Printer With Last 100 Cycle Storage Capacity.		
Dry Chamber	Use Vacuum Pump For Chamber Dry • Pre Vacuum • Post Vacuum • Chamber 100%		
Working Pressure	Variable 15 to 34 PSI		
For Over Pressure	Use Tra Fag Pressure Switch		
Working Temperature	121°C (up graded up to 150°C)		
Inner Chamber	Corrosion Resistant Stainless Steel SS 316L		
Outer Chamber	Corrosion Resistant 304 Grade Stainless Steel		
Finishing	Mirror Finished Inner And External Surface		
All Piping and Fittings	SS 304		
Door	Single or Double SS 304		
Insulation	Aluminum Silicate Which is the Best Insulating Material.		
Sterilization Period	25 to 30 Minutes		
For Water Level	Water Level Indicator		
When Water Level Low	Low Water Level Sensor Indicate in Display WATER LEVEL LOW.		
Steam Exhaust	5 to 7 Minutes		
Dry	5 - 10 Minutes		
Boiler	 Made of SS 304 Sheet Water Level Glass Tube x 2pcs Heating Element x 1pc Silicone Door Seal x 1pc 		
Chamber Condensate Line	Fitted with Steam Trap & NRV Valve		
Hydraulic Tested	On 66 psi		
Sensor	PT 100		
Safety Valve	For Excess Pressure Inside Jacket and Steam Generator.		
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes.		
Standard Accessories	 SS 304 Sterilizing Tray x 1pc Water Level Glass Tube x 2pcs Geating Element x 1pc Silicone Door Seal x 1pc 		

FEATURES & BENEFITS:

- Cylindrical Horizontal Autoclave
- Stainless Steel Interior And Exterior
- Triple Walled Construction
- Standard Temperature 121°c To 134°c
- Working Pressure 15 Psi To 34 PSI
- Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch
- Vacuum Breaker
 Steam Trap & NRV
- ISI Safety Valves
- Digital Temperature Controller With Timer And Alarm (buzzer)
- Automatic Water Level Cut-off
- Water Level Glass To Check Water Level
- SS Pressure & Compound Gauges
- ISIMarked Flange / Immersion Heating Elements



NEXX-002A

Silicone Door Seal x 1pc	Size of Cylindrical Horizontal Autoclave				
Fitted with Steam Trap & NRV Valve	INNER CHAMBER	INNER CHAMBER	CAPAC		
On 66 psi	IN MM	IN INCHES	ITY IN LTR	LOAD (KW)	
PT 100					
For Excess Pressure Inside Jacket and Steam Generator.	400 x 600 MM	16" X 24"	78	6-8	
Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes. • SS 304 Sterilizing Tray x 1pc • Water Level Glass Tube x 2pcs	400 X 900 MM	16" X 36"	120	8-10	
	400 X 1200 MM	16" X 48"	165	10-12	
	500 X 600 MM	20" X 24"	120	8	
Geating Element x 1pc Silicone Door Seal x 1pc	500 X 900 MM	20" X 36"	165	9	
220 Volts / 440 Volts (Three Phase)	500 X 1200 MM	20" X 48"	250	12	

NITIN KUMAR SHUKLA () +91-7007800569

Power Supply

Semi Automatic Rectangular Horizontal Autoclave Model Number: NEXX-003

NEXXGEN SURGICAL, Have Comprehensive Range Of Semi Automatic Rectangular Horizontal Autoclaves. Semi Automatic Rectangular Horizontal Autoclave Will Be Front Loading, With Microprocessor Based And In Horizontal Rectangular Shape With Triple Layer With Boiler. These High Pressure Front Loading Autoclaves Play Important Role In Sterilization Procedures Required In Hospitals, Pharmaceutical Companies, Biopharma And Microbiology Laboratories Etc.

This Autoclave Will Be Used In Sterilization Of Wrapped Instruments, Unwrapped Instruments, Glassware, Liquids, Metal Instruments, Culture Medium And Laboratory Purpose.

NEXXGEN SURGICAL Rectangular Horizontal Autoclave Units Are Made Of Stainless Steel To Ensure Corrosion Resistance Surface. Tipple Walled Construction Ensures Excellent Level Of Safety. Boiler Is Fitted Separately And Connected To The Chamber. Water Level Indicator Is Also Attached To The Boiler With Automatic Low Water Cut-off.

User Can Easily Set The Desired Temperature With The Help Of Digital Controller Fitted On The Front Side Of The Unit. We Provide Radial Locking Pressure Door With Silicon Gas Kit. We Use Automatic Pressure Control Switch And User Easly Can Set Temperature By Digital Controller. Our Machine Work With Pressure Between 15-34 Psi.

Our Horizontal Autoclaves Chamber And Boiler Are Hydraulic Tested On 66 PSI.

TECHNICAL SPECIFICATIONS

Model	NEXX-003		
Design	Rectangular Horizontal Type		
Construction	Triple Walled		
Automation Hardware	Microprocessor Based		
Display	3.6"(Inches). The Information Shows in Display Like Chamber Temperature, Chamber Pressure, Jacket Pressure and Timer Buzzer.		
Printer	Optional Feture (Thermal Printer)		
Dry Chamber	Use Vacuum Pump for Chamber Dry • Pre Vacuum • Post Vacuum • Chamber 100%		
Working Pressure	Variable 15 to 34 PSI		
For Over Pressure	Use TraFag Pressure Switch		
Working Temperature	121°C and 134°C (Up Graded Up to 150°C		
Inner Chamber	Corrosion Resistant Stainless Steel SS 316L		
Outer Chamber	Corrosion Resistant 304 Grade Stainless Steel		
Finishing	Mirror Finished Inner and External Surface		
All Piping and Fittings	SS 304		
Door	Single or Double Doors SS 304		
Insulation	Aluminum Silicate Which is the Best Insulating Material.		
Sterilization Period	25 to 30 Minutes		
For Water Level	Water Level Indicator		
When Water Level Low	Low Water Level Sensor Indicate in Display Water Level Low.		
Steam Exhaust	5 to 7 Minutes		
Dry	5 - 10 Minutes		
Boiler	 Made of SS 304 Sheet Water Inlet & Outlet Valves Water Level Glass Window Boiler Safety Valve With ISI Marked Water Immersion Type Industrial Heating Elements. 		
Chamber Condensate Line	Fitted With Steam Trap & NRV Valve		
Hydraulic Line	on 66 PSI		
Sensor	PT 100		
Safety Valve	For Excess Pressure Inside Jacket and Steam Generato		
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes		
Standard Accessories	 SS 304 Sterilizing Tray X 1Pc Water Level Glass Tube X 2Pc Heating Element X 1Pc Silicone Door Seal X 1Pc 		

FEATURES & BENEFITS:

- Rectangular Horizontal Autoclave
- Stainless Steel Interior And Exterior
- Triple Walled Construction Standard Temperature 121°c To 134°c
- Working Pressure 15 psi to 34 psi
 Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch Vacuum Breaker
- Steam Trap & NRV
 ISI Safety Valves
- Digital Temperature Controller With Timer And Alarm(buzzer)
- Automatic Water Level Cut-off
 Water Level Glass To Check Water Level
- SS Pressure & Compound Gauges
- ISI Marked Flange/immersion Heating Elements



Size of Rectangular Horizontal Autoclave

INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
450 x 450 x 900 mm	18"x18"x36"(1.5x1.5x3FT)	180	10KW
450 x 600 x 900 mm	18"x24"x36"(1.5x2x3FT)	250	12KW
600 x 600 x 900 mm	24"x24"x36"(2x2x3FT)	320	18KW
600 x 600 x 1200 mm	24"x24"x48"(2x2x4FT)	430	18KW
600 x 900 x 1500 mm	24"x36"x60"(2x3x5FT)	810	24-27KW
Customized Size			

NEXX-003

NITIN KUMAR SHUKLA () +91-7007800569

Power Supply

220 Volts / 440 Volts (Three Phase)

Fully Automatic Rectangular Horizontal Autoclave Model Number: NEXX-003B

NEXXGEN SURGICAL, have comprehensive range of Fully automatic Rectangular horizontal autoclaves. Fully automatic Rectangular horizontal autoclave will be front loading, with PLC/HMI Based and in Horizontal shape with triple layer with boiler. These high pressure front loading autoclaves play important role in sterilization procedures required in hospitals, pharmaceutical companies, biopharma and microbiology laboratories etc.

This Autoclave Will Be Used In Sterilization Of Wrapped Instruments, Unwrapped Instruments, Glassware, Liquids, Metal Instruments, Culture Medium And Laboratory Purpose.

NEXXGEN SURGICAL horizontal autoclave units are made of stainless steel to ensure corrosion resistance surface. Tipple walled construction ensures excellent level of safety. Boiler is fitted separately and connected to the chamber. Water level indicator is also attached to the boiler with automatic low water cut-off.

User can easily set the desired temperature with the help of digital controller fitted on the front side of the unit. we provide radial locking pressure door with silicon gas kit. we use automatic pressure control switch and user easly can set temperature by Digital Controller. our machine work with pressure between 15-34 PSI.

Our Horizontal Autoclaves Chamber And Boiler Are Hydraulic Tested On 66 PSI.

TECHNICAL SPECIFICATIONS

Model	NEXX-003B
Design	Rectangular Horizontal Type
Construction	Triple Walled
Automation Hardware	PLC / HMI
Display	7" (Inches) Touch Screen Milti-Colour. The Information Shows in Display Like Cycle Stage, Chamber Temperature, Chamber Pressure, Jacket Pressure, Along with the Information About Failures and Interrupts.
Printer	Thermal Printer With Last 100 Cycle Storage Capacity.
Dry Chamber	Use Vacuum Pump for Chamber Dry • Pre Vacuum • Post Vacuum • Chamber 100%
Working Pressure	Variable 15 to 34 PSI
For Over Pressure	Use TraFag Pressure Switch
Working Temperature	115°C and 138°C
Inner Chamber	Corrosion Resistant Stainless Steel SS 316L
Outer Chamber	Corrosion Resistant 304 Grade Stainless Steel
Finishing	Mirror Finished Inner and External Surface
All Piping and Fittings	SS 304
Door	Single or Double Doors SS 304 Hinges / Sliding
Insulation	Aluminum Silicate Which is the Best Insulating Material.
Sterilization Period	25 to 30 Minutes
For Water Level	Water Level Indicator
When Water Level Low	Low Water Level Sensor Indicate in Display Water Level Low.
Steam Exhaust	5 to 7 Minutes
Dry	5 - 10 Minutes
Boiler	 Made of SS 304 Sheet Water Inlet & Outlet Valves Water Level Glass Window Boiler Safety Valve With ISI Marked Water Immersion Type Industrial Heating Elements.
Chamber Condensate Line	Fitted With Steam Trap & NRV Valve
Hydraulic Line	on 66 PSI
Sensor	PT 100
Safety Valve	For Excess Pressure Inside Jacket and Steam Generator.
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes
Standard Accessories	 SS 304 Sterilizing Tray X 1Pc Water Level Glass Tube X 2Pc Heating Element X 1Pc Silicone Door Seal X 1Pc
Power Supply	220 Volts / 440 Volts (Three Phase)
eapp.,	

FEATURES & BENEFITS:

- Rectangular Horizontal Autoclave
- Stainless Steel Interior And Exterior
- Triple Walled Construction Standard Temperature 115°c To 138°c
- Working Pressure 15 psi to 34 psi
- Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch
- Vacuum Breaker Steam Trap & NRV ISI Safety Valves
- Digital Temperature Controller With Timer And Alarm(buzzer)
- Automatic Water Level Cut-off
- Water Level Glass To Check Water Level
- SS Pressure & Compound Gauges
- ISI Marked Flange/immersion Heating Elements



NEXX-003B

Size of Rectangular Horizontal Autoclave

INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
450 x 450 x 900 mm	18"x18"x36"(1.5x1.5x3FT)	180	12KW
450 x 600 x 900 mm	18"x24"x36"(1.5x2x3FT)	250	12KW
600 x 600 x 900 mm	24"x24"x36"(2x2x3FT)	320	18KW
600 x 600 x 1200 mm	24"x24"x48"(2x2x4FT)	430	18KW
600 x 900 x 1500 mm	24"x36"x60"(2x3x5FT)	810	24-27KW
Customized Size			

NITIN KUMAR SHUKLA (1) +91-7007800569

Sliding Door Fully Automatic Rectangular Horizontal Autoclave Model Number: NEXX-003BS

NEXXGEN SURGICAL, have comprehensive range of Fully automatic Sliding Door Rectangular horizontal autoclaves. Fully automatic Sliding Door Rectangular horizontal autoclave will be front loading and door automatically slid UP & DOWN, with PLC/HMI Based and in Horizontal shape with triple layer with boiler. These high pressure front loading Sliding Door autoclaves play important role in sterilization procedures required in hospitals, pharmaceutical companies, biopharma and microbiology laboratories etc.

This autoclave will be used in sterilization of wrapped instruments, unwrapped instruments, glassware, liquids, metal instruments, culture medium and laboratory purpose.

NEXXGEN SURGICAL horizontal autoclave units are made of stainless steel to ensure corrosion resistance surface. Tipple walled construction ensures excellent level of safety. Boiler is fitted separately and connected to the chamber. Water level indicator is also attached to the boiler with automatic low water cut-off.

User can easily set the desired temperature with the help of digital controller fitted on the front side of the unit. we provide radial locking pressure door with silicon gas kit. we use automatic pressure control switch and user easly can set temperature by Digital Controller, our machine work with pressure between 15-34 PSI.

Our Horizontal Autoclaves Chamber And Boiler Are Hydraulic Tested On 66 PSI.

TECHNICAL SPECIFICATIONS

Model	NEXX-003BS		
Design	Fully Automatic Sliding Door Rectangular Horizontal Type		
Construction	Triple Walled		
Automation Hardware	PLC / HMI		
Display	7" (Inches) Touch Screen Milti-Colour. The Information Shows in Display Like Cycle Stage, Chamber Temperature, Chamber Pressure, Jacket Pressure, Along with the Information About Failures and Interrupts.		
Printer	Thermal Printer With Last 100 Cycle Storage Capacity.		
Dry Chamber	Use Vacuum Pump for Chamber Dry • Pre Vacuum • Post Vacuum • Chamber 100%		
Working Pressure	Variable 15 to 34 PSI		
For Over Pressure	Use TraFag Pressure Switch		
Working Temperature	121°C (up Graded up to 150°C)		
Inner Chamber	Corrosion Resistant Stainless Steel SS 316L		
Outer Chamber	Corrosion Resistant 304 Grade Stainless Steel		
Finishing	Mirror Finished Inner and External Surface		
All Piping and Fittings	SS 304		
Door	Single or Double / Fully Automatic Sliding Door SS 304. (Door Slide Vertically UP & SOWN).		
Insulation	Aluminum Silicate Which is the Best Insulating Material.		
Sterilization Period	25 to 30 Minutes		
For Water Level	Water Level Indicator		
When Water Level Low	Low Water Level Sensor Indicate in Display Water Level Low.		
Steam Exhaust	5 to 7 Minutes		
Dry	5 - 10 Minutes		
Boiler	 Made of SS 304 Sheet Water Inlet & Outlet Valves Water Level Glass Window Boiler Safety Valve With ISI Marked Water Immersion Type Industrial Heating Elements. 		
Chamber Condensate Line	Fitted With Steam Trap & NRV Valve		
Hydraulic Line	on 66 PSI		
Sensor	PT 100		
Safety Valve	For Excess Pressure Inside Jacket and Steam Generator.		
Mechanical Time w / Alarm	Regulates Sterilization Time and Activates Buzzer When Sterilization Cycle Completes		
Standard Accessories	 SS 304 Sterilizing Tray X 1Pc Water Level Glass Tube X 2Pc Heating Element X 1Pc Silicone Door Seal X 1Pc 		
Power Supply	220 Volts / 440 Volts (Three Phase)		

FEATURES & BENEFITS:

- Fully Automatic Silding Door Rectangular Horizontal Autoclave
- Stainless Steel Interior And Exterior
 Triple Walled Construction
- Standard Temperature 121°c To 134°c
- Working Pressure 15 psi to 34 psi
- Silicon Jointless Rubber Door Gasket
- Automatic Pressure Control Switch
- Vacuum Breaker
 Steam Trap & NRV
 ISI Safety Valves
- Digital Temperature Controller With Timer And Alarm(buzzer)
- Automatic Water Level Cut-off
- Water Level Glass To Check Water Level
- SS Pressure & Compound Gauges
- ISI Marked Flange/immersion Heating Elements



Size of Rectangular Horizontal Autoclave

INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
450 x 450 x 900 mm	18"x18"x36"(1.5x1.5x3FT)	180	12KW
450 x 600 x 900 mm	18"x24"x36"(1.5x2x3FT)	250	12KW
600 x 600 x 900 mm	24"x24"x36"(2x2x3FT)	320	18KW
600 x 600 x 1200 mm	24"x24"x48"(2x2x4FT)	430	18KW
600 x 900 x 1500 mm	24"x36"x60"(2x3x5FT)	810	27KW
Customized Size			

NITIN KUMAR SHUKLA (1) +91-7007800569

Power Supply

220 Volts / 440 Volts (Three Phase)

Table Top Ethylene Oxide Gas Sterilizer Model Number : NEXX-004A

NEXXGEN SURGICAL, Table Top Ethylene Oxide Gas Sterilizers (ETO) are specifically designed to perform sterilization of temperature sensitive products. Our Sterilizercan be used both in Pharmaceutical and Parapharmaceutical industries and even in a wide range of medical and laboratory applications.

MAIN APPLICATIONS

- Syringes, Catheters, Cartridges for dialysis
- Bandages, Sutures
- Use by Laparoscopy Surgeon, Ophthal Surgeon
- Use by Ortho Surgeon, Pediatric Surgeon

• Plastic Articles

NEXXGEN SURGICAL manufacturing Ethylene Oxide Gas Sterilizers (ETO) as per the norms of the international standards and use branded parts and components. Material use in manufacturing only Stainless Steal(JINDAL) Mate and Mirror Finish both. Grade Use Only SS316 and SS304.

NEXXGEN Ethylene Oxide Gas Sterilizers Chamber tested under-600PSI Vacuum tested.

TECHNICAL SPECIFICATIONS

Model	NEXX-004A	
Design	Rectangular Horizontal Type	
Construction	Complete Stainless Steal	
Inner Chamber	Made by SS 316	
Outer and All Piping	Made by SS 304	
Automation Hardware	PLC / HMI	
Chamber Operating Pressure	Negative Pressure Cycle	
Chamber Testing Pressure	1.0 Kg / cm2	
Operating Temperature	Maximum: 65°C (chamber)	
Single Door	SS 304	
Door Locking	Manual / Automatic Pressure Locking System	
Gasket	Jointless Silicon Gasket	
Finishing	Mirror Finished Inner and External Surface Mate/Mirror.	
Insulation	Glass Wool Fiber 50mm	
Sterilization Period	6-7 Hours	
Heating	Heaters on all side of the Chamber for uniform heating throughout the load.	
Temperature Sensor	For accurate measuring and controlling of Temp. in the chamber	
Pressure Transmitter	For measuring and controlling the pressure-vacuum that the load is subjected to during removal and charging of gas.	
All Valves for Automation	JANATICS Pneumatic Solenoid Valves.	
Sensor	PT 100	
Heater	ESCORT	
Compressor	oil Free	
Vacuum Pump	Inbuilt	
	• SS 304 Loading Basket x 1 pc	
	Sealing Machine x 1 pcs	
Standard Accessories	• ETO Paper Roll 20cm x 1pc	
	• ETO Cartridges x 10pc	
	Air Blow Gun x 1pc	
Power Supply	220 Volts (single Phase)	

FEATURES & BENEFITS:

- Fully automatic One touch operated machine
- Thermal printer Inbuilt vacuum pump Inbuilt Air Compressor
- Rectangular chamber
- Stainless steel interior and exterior
- Standard temperature 35°C to 65°C
- Working pressure is -500 PSI.
- Silicon jointless rubber door gasket
- Manual/Automatic door lock
- Automatic Ethylene Cartridges puncher system.
- ISI marked flange/immersion heating elements

Small ETO Sterilizer (e.g 25-150Ltr) has been specially designed for Super Specialty and Multi-Specialty Hospitals capable for sterilizing all major medical devices.



Size of Rectangular Horizontal Autoclave

INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
310x310x310mm	12" x 12" x 12"	25LTR	1
310x310x600mm	12" x 12" x 24"	57LTR	2
310x310x900mm	12" x 12" x 36"	86LTR	3
310x310x1200mm	12" x 12" x 48"	115LTR	4
310x310x1350mm	12" x 12" x 54"	130LTR	4.5
310x310x1500mm	12" x 12" x 60"	145LTR	5

Ethylene Oxide Gas Sterilizer Model Number: NEXX-004B

NEXXGEN SURGICAL, Ethylene Oxide Gas Sterilizers (ETO) are specifically designed to perform sterilization of temperature sensitive products. Our Sterilizercan be used both in Pharmaceutical and Parapharmaceutical industries and even in a wide range of medical and laboratory applications.

MAIN APPLICATIONS

- Syringes, Catheters, Cartridges for dialysis
- Bandages, Sutures
- Use by Laparoscopy Surgeon, Ophthal Surgeon
- Use by Ortho Surgeon, Pediatric Surgeon

Plastic Articles

NEXXGEN SURGICAL manufacturing Ethylene Oxide Gas Sterilizers (ETO) as per the norms of the international standards and use branded parts and components. Material use in manufacturing only Stainless Steal(JINDAL) Mate and Mirror Finish both. Grade Use Only SS316 and SS304.

NEXXGEN Ethylene Oxide Gas Sterilizers Chamber tested under-600PSI Vacuum tested.

TECHNICAL SPECIFICATIONS

Model	NEXX-004B
Design	Rectangular Horizontal Type
Construction	Complete Stainless Steal
Inner Chamber	Made by SS 316
Outer and All Piping	Made by SS 304
Automation Hardware	PLC / HMI (4 & 7 INCHES)
Chamber Operating Pressure	Negative Pressure Cycle
Chamber Testing Pressure	1.0 Kg / cm2
Operating Temperature	Maximum: 65°C (chamber)
Single Door	SS 304
Door Locking	By 2Kg Pressure Automatic Locking System
Gasket	Jointless Silicon Gasket
Finishing	Mirror Finished Inner and External Surface Mate/Mirror.
Insulation	Glass Wool Fiber 50mm
Sterilization Period	6-7 Hours
Heating	Heaters on all side of the Chamber for uniform heating throughout the load.
Temperature Sensor	For accurate measuring and controlling of Temp. in the chamber
Pressure Transmitter	For measuring and controlling the pressure-vacuum that the load is subjected to during removal and charging of gas.
All Valves for Automation	JANATICS Pneumatic Solenoid Valves.
Sensor	PT 100
Heater	ESCORT
Compressor	Oil Free
Vacuum Pump	Inbuilt (OIL FREE)
Standard Accessories	 SS 304 Loading Basket x 1 pc Sealing Machine x 1 pcs ETO Paper Roll 20cm x 1pc ETO Cartridges x 10pc Air Blow Gun x 1pc
Power Supply	220 Volts (single Phase)

FEATURES & BENEFITS:

- Fully automatic One touch operated machine
- Thermal printer Inbuilt vacuum pump Inbuilt Air Compressor
- Rectangular chamber
- Stainless steel interior and exterior
- Standard temperature 35°C to 65°C Working pressure is -500 PSI.
- Silicon jointless rubber door gasket
 Manual/Automatic door lock

NEXX-004B

- Automatic Ethylene Cartridges puncher system.
- ISI marked flange/immersion heating elements

Small ETO Sterilizer (e.g 25-150Ltr) has been specially designed for Super Specialty and Multi-Specialty Hospitals capable for sterilizing all major medical devices.



Size of Rectangular Horizontal Autoclave

INNER CHAMBER IN MM	INNER CHAMBER IN INCHES	CAPAC ITY IN LTR	LOAD (KW)
310x310x310mm	12" x 12" x 12"	25LTR	1
310x310x600mm	12" x 12" x 24"	57LTR	2
310x310x900mm	12" x 12" x 36"	86LTR	3
310x310x1200mm	12" x 12" x 48"	115LTR	4
310x310x1350mm	12" x 12" x 54"	130LTR	4.5
310x310x1500mm	12" x 12" x 60"	145LTR	5

Industrial Ethylene Oxide Gas Sterilizer Model Number: NEXX-IETO-01

NEXXGEN SURGICAL, Industrial Ethylene Oxide Gas Sterilizers (Also Known as EO or ETO) are specifically designed to perform sterilization of temperature sensitive products Our Sterilizercan be used both in Pharmaceutical and Parapharmaceutical industries and even in a wide range of medical and laboratory applications.

MAIN APPLICATIONS

- Syringes, Catheters, Cartridges for dialysis
- Bandages, Sutures
- All disposable items used in hospitals

Plastic Articles

NEXXGEN SURGICAL manufacturing Industrial Ethylene Oxide Gas Sterilizers as per the norms of the international standards and use branded parts and components. Material use in manufacturing only Stainless Steal (JINDAL) Mate and Mirror Finis both. Grade Use Only SS316 and Ss304.

Industrial Ethylene Oxide Gas Sterilizers are two Type.

- First is manufactured in Triple walled with separate boiler system. Boiler System has used for uniform heat in chamber. Also provide dilution tank for ethylene Oxide.
- Second is manufactured in Double walled with Heater. Heater used for uniform heat in chamber.

NEXXGEN Ethylene Oxide Gas Sterilizers Chamber tested under-600PSI Vacuum tested.

FEATURES & BENEFITS:

- Fully automatic One touch operated machine
- Thermal printer
- Inbuilt vacuum pump
- Inbuilt Air Compressor
- Rectangular chamber
- Stainless steel interior and exterior
- Standard temperature 35°C to 65°C
- Working pressure is -300 to -400 PSI.

TECHNICAL SPECIFICATIONS

Model	NEXX-IETO-01
Design	Rectangular Horizontal Type
Construction	Complete Stainless Steal
Inner Chamber	Made by SS 316
Outer and All Piping	Made by SS 304
Automation Hardware	PLC / HMI
Chamber Operating Pressure	Negative Pressure Cycle
Chamber Testing Pressure	1.0 to 1.5 Kg / cm2
Chamber Operating Temp	Maximum: 65°C (Chamber)
Jacket Testing Pressure	3.0 Kg/cm2
Single / Double Door	SS 304
Door Locking	Manual Door Locking System
Gasket	Jointless Silicon Gasket
Finishing	Mirror Finished Inner and External Surface Mate/Mirror.
Insulation	Glass Wool Fiber 50mm
Sterilization Period	6-7 Hours
Heating	by Boiler
Temperature Sensor	For accurate measuring and controlling of Temp. in the chamber
Pressure Transmitter	For measuring and controlling the pressure-vacuum that the load is subjected to during removal and charging of gas.
All Valves for Automation	JANATICS Pneumatic Solenoid Valves.
Sensor	PT 100
Heater	ESCORT in Boiler
Water Pump	1HP Cromton
Vacuum Pump	Inbuilt
Standard Accessories	 SS 304 Loading Basket x 1 pc Sealing Machine x 1 pcs ETO Paper Roll 20cm x 1pc
Power Supply	220 Volts (single Phase)

- Silicon jointless rubber door gasket
- Manual door lock
- Ethylene Cylinder use for Ethylene Gas.
- ISI marked flange/immersion heating elements Use in Boiler



SIZE of ETO STERILIZER

	SIZE IN FEET W x H x L	CAPACITY IN LTR	VOLUMETRIC CAPACITY IN CFT	VACUUM PUMP	WATER CIRCULATION PUMP	ELECTRIC HEATERS LOAD
S.	2x2x4	430LTR	15	1HP	0.5HP	4KW
	2x2x5	550LTR	19	1HP	0.5HP	4KW
	3x3x6	1460LTR	51	1.5HP	0.5HP	6KW
	4x4x8	3500LTR	123	2HP	1HP	12KW
	5x5x10	6750LTR	238.5	3HP	1HP	12KW 18KW ^{to}
	6x6x16	15550LTR	549.5	5HP	2HP	18KW _{to} 24KW

BIO-MEDICAL WASTE TREATMENT CYLINDRICAL RECTANGULAR HORIZONTAL AUTOCLAVE Model Number: NEXX-BMW01

NEXXGEN SURGICAL, have comprehensive range of horizontal autoclaves for BIO-MEDICAL WASTE TREATMENT. Cylindrical/Rectangular Horizontal autoclave will be front loading and in Horizontal shape with triple layer with self steam generator or boiler. These high pressure front loading autoclaves play important role in BIO-MEDICAL WASTE TREATMENT procedures required in hospitals and in BIO-MEDICAL WASTE TREATMENT plants etc.

NEXXGEN SURGICAL BIO-MEDICAL WASTE TREATMENT autoclave units are made of stainless steel to ensure corrosion resistance surface. Tipple walled construction ensures excellent level of safety. Boiler is fitted separately and connected to the chamber. Water level indicator is also attached to the boiler with automatic low water cut-off.

User can easily set the desired temperature with the help of digital controller/HMI fitted on the front side of the unit. we provide radial locking pressure door with silicon gas kit. we use automatic pressure control switch and user easly can set temperature by Digital Controller, our machine work with pressure between 15-34 PSI.


Our BIO-MEDICAL WASTE TREATMENT autoclaves chamber and boiler are Hydraulic tested on 66 PSI.

TECHNICAL SPECIFICATIONS

Model	NEXX-BMW01	
Design	Cylindrical / Rectangular Horizontal Type	
Construction	Triple Walled	
Valve for Jacket to chamber Transfer, Chamber Release and for Post Vacuum or day mode	Use Heavy duty multiport or Valve	
Self Steam Generator System	SS 304 boiler with Heater	
Pre-Vacuum or High Vacuum	Yes	
Water in Boiler	Distilled/Normal Water	
Heater	Industrial Immersion Type water heater	
Working Pressure	Variable 15 to 34 PSI	
Working Temperature	121°C (up graded up to 135°C)	
Internal Chamber	Corrosion Resistant Stainless Steel SS 316	
External Chamber	Corrosion Resistant 304 Grade Stainless Steel	
Finishing	Mirror finished inner and external surface SS 304/316L	
All Piping	SS 316L	
Door	Single or Double Doors SS 304	
Insulation	Glass wool Fiber 100mm	
Sterilization Period	25 to 30 Minutes	
Steam Exhaust	5 to 7 Minutes	
Air removed during the vacuum cycle	Decontaminated by means of a HEPA filter (Class H13 or higher)	
Dry	5-10 Minutes	
Boiler	 Made of SS 304 sheet Water inlet & outlet valves Water level glass window Boiler safety valve With ISI marked water immersion type industrialheating elements. 	
Chamber Condensate Line	Fitted with steam trap & NRV valve	
Hydraulic Tested	on 66 psi	
Sensor	PT 100	
Temperature	Digital Display	
Temperature Indication	±1% accuracy	
Jacket Pressure	Use Pressure Gauge	
Chamber Pressure (for positive pressure)	Use Compound Gauge(for Positive & Negative Pressure)	
Pressure Indication	±1.6% accuracy	
Safety Valve	For excess pressure inside Chamber, jacket and steam generator.	
Mechanical Time w / Alarm	Regulates sterilization time and activates buzzer when sterilization cycle completes	
Standard Accessories	 SS 304 Sterilizing tray x 1pc Water level glass tube x 2pcs Heating element x 1pc Silicone door seal x 1pc 	
Optional Accessories	PLC operation Vacuum drying Printer Auto add water Door interlock Automatic Pressure control switch water cut-off Temperature indicator Caster wheels	
Power Supply	220 Volts/440 Volts (Three Phase)	

FEATURES & BENEFITS:

- Cylindrical/Rectangular Horizontal Autoclave
- Stainless steel interior and exterior
- Triple walled construction
- Self-locking safety door cannot be opened when chamber is under pressure
- Standard temperature 121°C to 134°C
- Working pressure 15 psi to 34 psi
- Silicon jointless rubber door gasket
- Automatic pressure control switch
- Vacuum breaker Steam trap & NRV ISI safety valves
- Digital Temperature controller with Timer and alarm (buzzer)
- Automatic water level cut-off
- Water level Glass to Check Water Level
- SS Pressure & Compound gauges
- ISI marked flange / immersion heating elements



Size of Cylindrical/Rectangular Horizontal Autoclave INNER CHAMBER IN MM INNER CHAMBER IN INCHES CAPACITY IN LTR LOAD(KW) 400 X 600 MM 16" X 24" 78 6-8 400 X 900 MM 16" X 36" 120 8-10 400 X 1200 MM 16" X 48" 165 9 500 X 600 MM 20" X 24" 120 8 500 X 900 MM 20" X 36" 165 9 500 X 1200 MM 20" X 48" 250 12 22" X 52" 550 X 1300 MM 310 18 INNER CHAMBER IN MM INNER CHAMBER IN INCHES CAPACITY IN LTR LOAD(KW) 450 X 450 X 900 MM 18" X 18" X 36" (1.5 X 1.5 X 3FT) 180
 450 X 600 X 900 MM
 18" X 24" X 36" (1.5 X 2 X 3FT)

 600 X 600 X 900 MM
 24" X 24" X 36" (2 X 2 X 3FT)
 250 320 600 X 600 X 1200 MM 24" X 24" X 48" (2 X 2 X 4FT) 450 600 X 900 X 1500 MM 24" X 36" X 60" (2 X 3 X 5FT)

10KW

12KW

18KW

18KW

24-27KW

810

Customized Size.....

Mortuary Chamber Model Number: NEXX-MC01

NEXXGEN SURGICAL, is manufacturer of Mortuary Chamber for storing cadaverous under low temperature cool conditions in order to prevent decomposition. Mortuary chambers have multiple names such as dead body freezer, mortuary cabinets, mortuary freezers, dead body storage freezer and body storage freezer. Mortuary Chamber Cabinet are double walled with outer made of SS 304 and inner walls made of SS 316. Tray is made of SS 304 to carry dead bodies. Trays of our Dead Body Cabinet travel on rollers fitted with ball bearings for easy movement. We provide hermetically sealed compressor Serves Refrigeration Unit and air conditioner unit use VOLTAS, BLUE STAR etc.

We provide digital microprocessor temperature indicator cum Controller to show temperature. Temperature range of morturay chmaber are-2°C to 10°C AND 8°C to-10°C CUSTOMISED.

Mortuary cabinets have a wide application and are used in hospitals, airports, research centers and many other agencies that may have a requirement to preserve a human body before being used in research or transported to destination.

The primary use of a mortuary chamber is to store and prevent a body from decomposing rapidly. Sometimes, it happens that a body has to be kept for long duration because of unexpected circumstances.

Technical specification of 4 Body Mortuary Cooling Chamber:

- Chamber storage temperature
- (-2°L to 10°C AND 8°C to -10°C Cooling)
- Inner/outer cabinet is made of stainless steel grade. ss304. (Sturcly construction)
- Puff insulation 100mm
- Microprocessor based temperature controller
- The cooling chamber system should be fans shut off & Interior proper LED lighting ON during door opening.
- Audio visual alarm for high & low temperature.
- Dead bodies trays with handles and travel roller fitted with ball bearing.
- The floor constructed out hardened and rain forced steel sheets containing spills for easy cleaning & hygiene (proper drain provide)
- Energy efficient & low maintenance

TECHNICAL SPECIFICATIONS

Model	NEXX-MC01
Design	Rectangular Horizontal Type
Construction	Double Walled
Working temperaturer	-2°L to 10°C AND 8°C to -10°C
Internal Chamber	Stainless Steel SS 316
External Chamber	Stainless Steel 304 Grade
Tray	Stainless Steel 304 Grade
Finishing	Mirror finished inner and external surface
Door Custruction	SS 304
Insulation	Puf 60mm
Sterilization Period	25 to 30 Minutes
Temperature Sensor	PT 100
Display	LCD/LED
Refrigerant (Condenser)	Cfc free
Air Circulation (Condenser)	Efficient condenser with automatic evaporating system (condensate) Forced
Controller (temperature)	PID
Body Capacity	1/2/3/4/6 (Customised)
Power Supply	220 V & 50/60 Hz
Cluster Wheel	Lockable

- Front panel of the units comprises of on/off switches coolingand main indication lamps, temperature controller, voltmeter, ammeter
- Interior led lighting Casters Swivel licking castors
- Superior engineered frame for maximum rigidity and low
- noise level <60dB
- Strong and reliable with powerful quick response quick response refrigeration, and exclusive ling life process is an
- unique feature of mortuary cabinet
- Heavy duty air cooled non CFC compressor Auto defrosting system
- ALL control like 1 pc, hpc, volt guard protection etc should be provide

The condensing unite & outdoor unit is around 10-20 feet

· Warranty period minimum one year



NEXX-MC01

Laminar Air Flow Model Number: NEXX-LAF01

NEXXGEN SURGICAL, are leading manufacturer, supplier and exporter of laminar air flow. we manufacturing Laminar Air Flow such as Vertical Laminar Air Flow and Horizontal Laminar Air Flow. A laminar air flow workstation is a closed cabinet fitted with HEPA filtered air flow system.

Laminar means unidirectional constant flow of air with almost no or minimal turbulence. The air flow velocity remains between 0.3 m/s to 0.5 m/s. The purpose of using such workstations in laboratory is to create particle and bacteria free working environment to carry out specialized work. As these units discharge air towards user, they provide no personal protection but product protection from room contaminants

Laminar Flow Cabinets, or laminar air flow cabinets as they are also known, are normally made of stainless steel with no gaps or joints thereby preventing the build-up of bacteria from collecting anywhere in the working zone.

Laminar Flow Cabinets are also known as clean benches because the air for the working environment is thoroughly cleaned by the precision filtration process.

HORIZONTAL LAMINAR AIR FLOW Cabinets receive their name due to the direction of air flow which comes from above but then changes direction and is processed across the work in a horizontal direction. The constant flow of filtered air provides material and product protection.

VERTICAL LAMINAR AIR FLOW Cabinets function equally well as horizontal Laminar Flow Cabinets with the laminar air directed vertically downwards onto the working area. The air can leave the working area via holes in flow cabinets can provide greater operator protection. the base. Vertical


HEPA FILTER filtered air through pre filter is forwarded towards HEPA FILTER. HEPA FILTER offers 99.99% efficiency for 0.3 micron particulates.

LAMINAR AIR FLOW APPLICATIONS:

- Life science research
- Microbiology
- Pathology lab
- Tissue culture

TECHNICAL SPECIFICATIONS

Model	NEXX-LAF01
Design	Rectangular Horizontal Type
Construction	SS 304 / Powder coated MS / GI sheet
Air flow direction	Horizontal/ Vertical
Pre-Filter	10 microns (washable)
HEPA filter	0.3 Microns(not washable it will change in 12 to 18 months)
Air velocity	0.45 m/s to 0.65 m/s
Front Door	Manual made by Acrylic transparent seat
Finishing	Mirror finished inner and external surface
Door	Single Door
Power Supply	220 volts / 50 Hz



- Musroom Spam Lab
- Laboratory
- Electronic and industrial sectors



Size of Laminar Air Flow

NEXX-LAF01

SIZE IN MM W X H X D	SIZE IN FEET W X H X D	CAPAC ITY IN LTR	HEPA FILTER SIZE/NO.	NO. OF PRE- FILTER		
600X600X600	24" X 24" X 24"	216LTR	2FT X 2FT X 6Inch / #1	#1		
900X600X600	36" X 24" X 24"	325LTR	3FT X 2FT X 6Inch / #1	#2		
1200X600X600	48" X 24" X 24"	430LTR	4FT X 2FT X 6Inch / #1	#2		
1800X600X600	72" X 24" X 24"	650LTR	3FT X 2FT X 6Inch / #2	#3		
2400X600X600	96" X 24" X 24"	870LTR	4FT X 2FT X 6Inch / #2	#4		

Hot Air Oven Model Number: NEXX-HAO01

NEXXGEN SURGICAL, are leading manufacturer, supplier and exporter of Hot Air Oven. The construction of hot air oven is double walled and made of 304 grades stainless steel with insulation of PUF for efficient thermal loss. In economical models, exterior chamber is made of powder coated MS, while 304 stainless steel is used as an option. Inner chamber is completely made of stainless steel 304, which is easy to clean and maintain. Shelves are made of painted steel and can be easily adjusted inside the chamber with ribs created inside the interior chamber walls.

Hot Air Oven temperature range 5°C above ambient to 250°C with accuracy of ± 2 °C. Temperature is controlled by digital microprocessor based PID temperature controller. For uniform rotation of heat inside the chamber noise-less fan is used ased for air circulation. Each unit is fitted with over temperature protection device (safety thermostat), which keeps your oven and samples safe in case of controller failure.

Hot Air Oven doors are made in double walled with PUF Insullation and Hinge Type.

VENTILATION For safe removal of hot gases and fumes air ventilators are provided near the top of the both sides

LAMINAR AIR FLOW APPLICATIONS:

- It uses dry heat for sterilizing
- Has a digital thermostat to regulate temperature
- Double wall design to retain heat and reduce energy output
- No need for extra water during the process

TECHNICAL SPECIFICATIONS

Model	NEXX-HSO01
Model	
Design	Horizontal Type
Construction	Double Walled
Working Temperature	5°C to 250°C
Temperature Accuracy	+/-2°C
Temperature Uniformity	+/-1°C
Internal Chamber	Corrosion Resistant Stainless Steel SS 304
Finishing	Mirror finished inner and external surface
Outer Body	SS 304/ MS
All Trays	SS 304
Door	Single Doors double walled
Insulation	PUF 60MM
Controls	PID Controller
Temp Display	LED Display
Sensor	PT 100
Heating Element	ISI marked industrial heating elements.
Safety device	Over temperature protectionElectric leakage breaker
Power Supply	220 Volts

- Smaller in size than the autoclaves
- Highly efficient and use less heat
- Used in laboratories to sterilize operating instruments
- Use in Mushroom Lab
- Use in Hospitals



SIZE OF HOT AIR OVEN

NEXX-HSO01

SIZE IN MM W X H X D	SIZE IN FEET W X H X D	CAPAC ITY IN LTR	LOAD (KW)	SHELVE S(NO's)		
300 x 300 x 300 mm	12" x 12" x 12"	30	1.5	2		
355 x 355 x 355 mm	14" x 14" x 14"	45	1.5	2		
455 x 455 x 455 mm	18" x 18" x 18"	95	2	2		
455 x 455 x 605 mm	12" x 12" x 24"	125	2.2	3		
605 x 605 x 605 mm	24" x 24" x 24"	225	2.5	3		
605 x 455 x 910 mm	12" x 18" x 36"	252	2.5	3		
605 x 605 x 910 mm	24" x 24" x 36"	336	3	3		

Bacterial Incubator Model Number: NEXX-BOD01

NEXXGEN SURGICAL, are leading manufacturer, supplier and exporter of Bacteriological Incubator. The construction of Bacteriological Incubator is double walled and made of 304 grades stainless steel with insulation of PUF for efficient thermal loss. Double door. Outer metalic double walled and inner glass door. In economical models, exterior chamber is made of powder coated MS, while 304 stainless steel is used as an option. Inner chamber is completely made of stainless steel 304, which is easy to clean and maintain. Shelves are made of SS 304 and can be easily adjusted inside the chamber with ribs created inside the interior chamber walls.

Bacteriological Incubator heating elements are designed with law wattage density to work on black heat. The heaters are in coil form and supported with refractory, covers complete bottom to give uniform heating. Specially designed thermostat gives better temperature accuracy. A Digital Temperature Controller can be used for temperature control.

The temperature range of our bacterial incubators is ambient+10°C to 60°C, which is controlled through microprocessor based PID temperature controller and PT100 RTD sensor. This controller displays both Set Value (SV) and Process Value (PV). PSI.

Ideal temperature for incubation is 37°C

LAMINAR AIR FLOW APPLICATIONS:

- Has a digital thermostat to regulate temperature
- Double wall design to retain heat and reduce energy output
- No need for extra water during the process

- Smaller in size than the autoclaves
- Highly efficient and use less heat
- Used in laboratories to sterilize operating instruments
- Use in Mushroom Lab
- Use in Hospitals

TECHNICAL SPECIFICATIONS

Model	NEXX-BOD01
Design	Square/ Rectangular Horizontal Type
Construction	Double Walled
Working Temperature	Ambient + 5°C to 60°C
Temperature Accuracy	0.5% °C
Temperature Uniformity	1% of set value
Internal Chamber	Corrosion Resistant Stainless Steel SS 304
Finishing	Mirror finished inner and external surface
Outer Body	SS 304/ MS
All Trays	SS 304
Door	Single Doors double walled
Insulation	PUF 60MM/ Glass Wool 65 mm.
Controls	Digital Display Electronic Controller.
Ventilator	Adjustable type (Aluminum)
Sensor	PT 100
Heating Element	ISI marked industrial heating elements.
Safety device	 Thermostat or Blind Controller ON/OFF cycle. Microprocessor based PID Controller Forced air circulation
Power Supply	220 Volts



SIZE OF BACTERIOLOGICAL INCUBATOR

SIZE IN MM W X H X D	SIZE IN FEET W X H X D	CAPAC ITY IN LTR	SHELVE S(NO's)
300 x 300 x 300 mm	12" x 12" x 12"	30	2
355 x 355 x 355 mm	14" x 14" x 14"	45	2
455 x 455 x 455 mm	18" x 18" x 18"	95	2
455 x 455 x 605 mm	12" x 12" x 24"	125	3
605 x 605 x 605 mm	24" x 24" x 24"	225	3
605 x 455 x 910 mm	12" x 18" x 36"	252	3
605 x 605 x 910 mm	24" x 24" x 36"	336	3

Ultrasonic Cleaner Model Number: NEXX-UC01

NEXXGEN SURGICAL, is recognized as one of the well-established manufacturers and exporters of ultrasonic cleaners for hospitals & laboratories. Ultrasonic cleaner is microprocessor based, efficient and effective equipment widely used for cleaning purposes in hospitals. The machine are manufactured by using high grade materials and are widely appreciated for their features like user friendly operations, excellent endurance capacity and optimum cleaning quality.

There are in all 3 cycles of operation in the ultrasonic tank:

- Water In, in which the water starts flowing in the tank for a preset time
- Ultrasonic, in which the objects are cleaned by ultrasonic cavitations
- Ultrasonic with Swinging, in which the basket containing the objects starts swinging in the ultrasonic tank and are thoroughly cleaned by dual action of ultrasonic cavitations and swinging movements.



TECHNICAL SPECIFICATIONS

Model	NEXX-UC01
Design	Rectangular Horizontal Type
Working Chamber Material	SS 316
Outer Material	SS 304
Automation Grade	Semi / Automatic
Frequency	30+/-3 Khz
Thermostate	Digital
Finishing	Mirror finished inner and external surface
Heater	1000 WATTS, CARTRIDGE TYPE
Average Power Output	500 WATTS



30KG COMMERCIAL WASHING MACHINE











60KG COMMERCIAL WASHING MACHINE

Description

PRODUCT DETAILS:

Capacity Up to	: 60 kg
Area in square-mm	: 910x1220
Basket Lenght(mm) :1220
Basket width(mm)	: 910
Frequency	: 50-60Hz
Power Source	: Electric
Country of Origin	: Made in India
Warranty	: 1 years
Convenience	:to do laundry together
Efficiency	; 60 kg
Automation Grade	: Semi-Automatic
Phase	: 3
Voltage	: 410
Motor Power	: 3 hp
l Deal In	: New Only
Loading Type side	: side Loading
Usage/Application	: Related to laundry
Basket size in inch	: 36/48

LANGUAGE

SURGICALS

- English
- Hindi


SKILLS SUMMARY

- Low water consumption
- Low Steam consumption
- Less time in Process
- More durability
- More Loading capacity

NITIN KUMAR SHUKLA () +91-7007800569

PERFORMANCE

Motor Brand	: N/A
Volt	: 410V
Rpm	: 32
VFD(Optiona	: Adjustable Drum Rotation 0-32 RPM
Forward/reverse	: Adjustable 0-60 Seconds
Beaters	:To improve wash Quality
Washing Time	: 0-60 min controlled by Timer
Temperature	: N/A
Control panel	: Manually Operated by Switch
Water level	: N/A
	Volt Rpm VFD(Optiona Forward/reverse Beaters Washing Time Temperature Control panel





CATALOG FOR

COMMERCIAL WASHING MACHINE

- MODEL NO. PSV100, PSV60 PSV50, PSV40, PSV30,
- Commercial washing machines play an important role in washing laundry clothes and washing one or more clothes.
- semi automatic washing machine
- Semi automatic comes at a higher price than automatic machine.
- Commercial washing machine drums are made from 1.2 mm (18 gauge) thick steel to withstand wear and tear over time.
 8mm iron is applied on the side plate of the commercial washing machine drum, so that the drum does not tear quickly.
- The biggest problem in commercial washing machines is that their metal melts, keeping this in mind, our company gives you a steel in the lower body in this price range so that it does not melt quickly.



- More Loading capacity
- Low water consumption
- Low Steam consumption
- Less time in Process
- More durability

	Basket width(mm)	Basket Lenght(mm)	square-	basket size in inch	Metal Type	Motor(HP)	MOTOR RPM	Panel Type	Power Source	Frequency
100	1020	1525	1020x1525	40/60	Stainless Steel	3	960	Semi-Automatic	Electric	50-60Hz
60	910	1220	910x1220	36/48	Stainless Steel	3	1440	Semi-Automatic	Electric	50-60Hz
50	910	1120	910x1120	36/45	Stainless Steel	3	1440	Semi-Automatic	Electric	50-60Hz
40	810	1110	810x1110	32/45	Stainless Steel	2	1440	Semi-Automatic	Electric	50-60Hz
30	660	910	660x910	26/36	Stainless Steel	1	1440	Semi-Automatic	Electric	50-60Hz
18	560	660	560x660	22/26	Stainless Steel	1	1440	Semi-Automatic	Electric	50-60Hz

Optional:

- (a)Semi automatic Based Controller
- (b) able to wash all types of clothes
- (c) Manual and automatic control panel
- Panel-Semi automatic
- Body Matel-S.S. J4 (202/304)
- Durm matel-S.S. J4 (202/304)
- capacity KG 30,40,50,60,100







CATALOG FOR

TUMBLE DRYER MACHINE







FEATURES

Capacity (KG)	Basket width(mm)	Basket Lenght(mm)	Area in square-mm	Metal Type	Motor (HP)	RPM	Power Source	Frequency
60	920	1260	920x1260	Stainles s Steel	2	1440	Electric	50-60Hz
50	910	1220	910x1220	Stainles s Steel	2	1440	Electric	50-60Hz
30	800	910	800x910	Stainles s Steel	2	1440	Electric	50-60Hz
25	700	910	700x910	Stainles s Steel	2	1440	Electric	50-60Hz

LAUNDRY RELATED MACHINE



15 INCH HYDRO DRYER MACHINE

Description

PRODUCT DETAILS:

PERFORMANCE

Motor Brand

Temperature

Volt

Beaters

: N/A

: 410V

: N/A

:To improve wash Quality

Capacity Up to	: 12 kg
Area in square-mm	: 300x380
Basket Lenght(mm):380
Basket width(mm)	: 380
Frequency	: 50-60Hz
Power Source	: Electric
Country of Origin	: Made in India
Warranty	: 1 years
Convenience	:to do laundry together
Efficiency	; 12 kg
Automation Grade	: Direct
Phase	: 3
Voltage	: 410
Motor Power	: 0.5 hp
l Deal In	: New Only
Loading Type side	: Top Loading
Usage/Application	: Related to laundry
Basket size in inch	: 12/15

<u>LANGUAGE</u>

- English
- Hindi

Ų


4

J

Ο

SKILLS SUMMARY

- •
- •
- Less time in Process
- More durability
- More Loading capacity







24 INCH HYDRO DRYER MACHINE

18 INCH HYDRO DRYER MACHINE








































NITIN KUMAR SHUKLA Proprietor

+91-7007800569



🔯 nitinshukla.shukla2@gmail.com