



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY  
 905, Panorama complex, No. 104001, Sector 100, Gurgaon, Haryana,  
 Aspians road, Sector 100, Haryana, India  
 P.T. KRUPA MEDICALS & SERVICES LTD., 17, Bawali Road,  
 E-mail: [info@krupamedical.com](mailto:info@krupamedical.com)

## Cavitation + Tripolar RF + Monopolar RF + Viscuatt™



### Features

- Patients and non-surgical way to remove fat and cellulite.
- Reduces fat in thighs, hips, abdomen, back, arms, chin etc.
- Tighten and tone skin.
- Body reshaping.
- Enhancing muscle flexibility.

### Technical Specification

Power supply input: 110V/220V 50Hz/60Hz

Power: 100W

40K cavitation head

Power supply output: 110V

Frequency: 40KHz

Power: 25W

Vacuum+Bipolar RF head:

Power supply output: 110V

Frequency: 4500Hz

Power: 30W

Tripolar RF head:

Power supply output: 110V

Frequency: 4500Hz

Power: 30W

Face RF head:

Power supply output: 110V

Frequency: 4500Hz

Power: 30W

### Slimming for below body parts



### Accessories

Monopolar RF head/Tripolar RF) x 1

40KHz Cavitation Head x 1

Vacuum Bipolar Head x 1

Bipolar RF Head for face x 2

Power cord x 1

Petit Bag x 1 Small x 2

Certified By:



ISO 13485 (Medical Standard)


**KRUPA**

AN ISO 9001:2008 CERTIFIED COMPANY

305,Panorama complex, Nr Navjivan Press, B/h Gujrat vidhyapith,  
Ashram road, Income tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com



**Cavitation + Tri Polar Rf  
+  
Mono Polar Rf + Vacuum**



Inch Loss



Tummy Tuck



Body Shaping

#### Specifications

TFT LCD Display	8 inch color touch screen
Technology	Cavitation+TriPolar RF+Monopolar RF+Vacuum
Cavitation Frequency	40KHz , 0 to 50W standard. 0 to 100W for option
TriPolar RF,Monopolar RF	1MHz, 10 MHz (output frequency) ; 100J (Maximum output power) 1MHz, 10 MHz(output frequency) ; 100J (Maximum output power)
Vacuum Pressure	1000Pa
Handpieces	1). Lipolysis RF head 2) Tripolar RF head (one big, one small ) 3) Ultrasonic cavitation explode head 4) Vacuum liposuction head
Voltage	220v AC/110 AC
Pack Size	65cmx55cmx55cm (LxWxH)
Net Weight	25kg
Gross Weight	34kg



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

KRUPA

305, Panorama complex, Nr. Navjivan Press, B/h Gujarat vidhyapith,  
Ashram road, Income tax, Ahmedabad-380014.  
Ph-9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Cosmosurge

### 4.1MHz Ablative-Nonablative With Hair Removal



#### Applications

- ⇒ Compound Nevus
- ⇒ Epidermal Nevus
- ⇒ Depressed Scars
- ⇒ Warts,Mole,Tags
- ⇒ Actenickeratosis
- ⇒ Granulom
- ⇒ Birthmark

#### Specification

- + Power Required 220-240v/50Hz
- + Insulated
- + 4.1 Mhz
- + 70%cut 30%coag.
- + 60%coag 40%cut.
- + Available.
- + Digital

#### Features

- ↗ 3 in 1 Ablative & Non-Ablative RF with Hair Removal
- ↗ Elegant table top model
- ↗ Patient isolation transformer for complete safety
- ↗ No need of ground contact or skin contact of patient plate
- ↗ Pain less hair removal probe with conductivity gel
- ↗ No need anesthesia for Warts,Mole or Tags





# KRUPA MEDI SCAN

AN ISO 9001-2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Professional Microdermabrasion

### High performance 4 in 1 Technology



Crystal



+ Diamond



+ Meso



+ Oxygen Spray

#### Technical Specification



- Hand piece material : As per medical grade
- Hand Piece : Pencil type
- Tips material : Standard medical grade(can be sterilized)
- Diamond probe : Pencil type
- Diamond tips : 3 Pcs as per medical standard
- Dehumidifier : with Hydro beads
- Delivery angle : 90 degree
- Exfoliation agent : Pure aluminium oxide sterilized crystals
- Particle distribution : 100/140/180 Microns
- Containers : Fresh supply (3) jars
- Vacuum pump : Waste collection jar (1)
- Powder flow : 2 nos imported noise free pump
- Vacuum : Automatic as per vacuum adjustment
- Operating system : 26 inch above
- Electrical supply : Micro processor base
- Meso therapy : S.M.P.S 220V/50Hz (In build spike guard protection)
- Meso therapy frequency : Transderm technology (Ultrasound + Iontophoresis)
- Meso therapy current : 1MHz
- Oxygen air output : 1mA to 30 mA
- Air pressure : 8LPM approx
- Dimension (cm) : 2.5bar approx
- Weight : 36 L x 30 W x 26 H
- Weight : 20kg approx



KRUPA

# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305,panoramic complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## *Derma Care*

### *Applications*

- Active acne scar revision
- Super facial pigmentation
- Un even skin tone
- General dermabrasion
- Face and body peeling lined
- Wrinkled or sun damaged skin
- Age spot from the face
- Remove black & white heads
- Skin Polishing.



### *Features*

- Virtually no side effects
- Effective on all skin color & types
- No anesthesia or recovery time required
- Immediate result & continuing improvement
- Patient return immediately to daily activities
- can be performed quickly (about 15 to 20 Minutes)
- 2 in 1 model - Crystal+Diamond



Before



After



### *Technical Features*

- Noise less
- Imported vacuum pump
- Vacuum capacity: 25 inch above
- Power- AC 220/240v-50Hz
- Crystal flow: Manual as per requirement
- Delivery Angle: 90 degree
- Media Material: Pure aluminium oxide sterilized crystals
- Crystal probe tip: as per medical grade
- Hand piece material: medical standard
- Particle distribution: 80% + 100 microns
- Weight: 9.5 kg approx
- Size 32x31x21cm

Certified By :



ISO 13485 (Medical Standard)



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Derma Chair



**Eco**

**Delux**

**Super Delux**

A specially designed export quality derma chair with high quality cushion & superior furnishing for comfort & fully flexible recliner chair 180 degree rotation convenient for doctor

### Features

- Luxurious styling multi position section couch suitable for all dermatology,Cosmetology procedures
- Imported DC German motor ensures smooth & jerk free height adjustment
- Fully motorize chair operated by foot switch or hand control remote for UP,DOWN,Front & BACK movements
- Instead of moving back side,doctor may recline the chair back on his front & check the patient very easily
- The centrally located 180degree rotation control maximizes space saving
- Easy to operate & practically maintenance free

	Eco	Delux	Super Delux	Dimension	Eco	Delux	Super Delux
Power supply	230v,50Hz,AC	230v,50Hz,AC	230v,50Hz,AC	Total length after recline	60.0"	67.0"	72.0"
Motor type	Imported Motor	Imported DC German Motor	Imported DC German Motor	Total length normal position	26.0"	28.0"	30.0"
Operating control	Foot switch or Hand control	Foot switch or Hand control	Foot switch or Hand control	Total length b/w two handle	22.0"	24.0"	26.0"
Chair rotation	180 degree	180 degree	180 degree	Total height	46.0"	50.0"	54.0"
Vertical traveling	200mm	280mm	280mm				
Power consumption	100 w	100 w	100 w				
Operating mode	Front,Back,Up,Down,Rotation	Front,Back,Up,Down,Rotation	Front,Back,Up,Down,Rotation				
Weight capacity	100kg Approx	120kg Approx	150kg Approx				



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305,panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## **Derma Gold**



3 Different Outlet  
100/140/180 Microns

### Applications

- Active acne scar revision
- Hyperpigmentation
- Stretch marks
- Dull,Congested, Thick or Blemished skins
- General Dermabrasion
- Face & body peeling lined
- Wrinkled or sundamaged skin
- Tattoo removal
- Age spot from the face

### Features

- ☒ Virtually no side effects
- ☒ Effective on all skin color & types
- ☒ No anesthesia or recovery time required.
- ☒ Immediate result & continuing improvement
- ☒ Patient return immediately to daily activities
- ☒ can be performed quickly (about 15 to 20 Minutes)
- ☒ 2 in 1 model- Crystal+Diamond
- ☒ 3 different outlets 100/140/180 microns
- ☒ Non blocking systems
- ☒ Automatic crystal flow

### Result after 7 sittings



Before

After

### Technical Features

- |                         |  |
|-------------------------|--|
| ◆ Hand piece material   | : As per medical grade                       |
| ◆ Hand piece            | : Pencil type                                |
| ◆ Tips material         | : Standard medical grade (can be sterilized) |
| ◆ Diamond probe         | : Pencil type                                |
| ◆ Diamond tips          | : 3 Pcs as per medical standard              |
| ◆ Delivery angle        | : 90 degree                                  |
| ◆ Exfoliation agent     | : Pure aluminium oxide sterilized crystals   |
| ◆ Particle distribution | : 100 / 140 / 180 Microns                    |
| ◆ Containers            | : Fresh supply (3) jars                      |
| ◆ Vacuum pump           | : Waste collection (1) jar                   |
| ◆ Powder flow           | : 1Nos imported noise free pump              |
| ◆ Vacuum                | : Automatic as per vacuum adjustment         |
| ◆ Operating system      | : 25 inch                                    |
| ◆ Electrical supply     | : Micro processor base                       |
| ◆ Dimension (cm)        | : SMPS 220V 50Hz (In build spiker)           |
| ◆ Weight                | : L: 36 x W: 30 X H: 26                      |
|                         | : 12kg approx                                |



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kmts@krupamediscan.com

## Automatic Micro needling Therapy Systems (AMTS)

### DERMA PEN



Auto-MTS Therapy is Derma Rejuvenation Therapy Auto-MTS is Non-surgical skin Rejuvenation treatment which utilizes Micro Rolling Needle.

Change the future of derma rejuvenation Derma Pen is a solution for the meso-roller

- => portable Wireless
- => Automatic micro needle
- => Therapy care system
- => low pain
- => depth control available
- => competitive price
- => rechargeable



#### Application

- => Skin rejuvenation
- => Enhancement of elasticity
- => Wrinkle reduction
- => Acne reduction
- => Stretch mark reduction
- => Stimulation of hair growth
- => Whitening
- => Administration of active ingredients

#### Features

- => A system with multi-acupuncture needles to stimulate the skin.
- => Effective use in hospitals as AMTS(Automatic Micro-needle Therapy System)
- => Convenient use with rechargeable batteries
- => The whole family can use it with a simple needle replacement
- => Safety and effectiveness as a medical device approved by korea
- => Certified by KFDA, CE, FDA, 3patents
- => 3 patents received for excellence in technology

#### Basic Specification

- => Multi-needling system stimulating skin
- => Depth adjustable (0.25mm~2.0mm) only with 1 needle
- => Much decreased pain compared to meso roller
- => Excellent improving effect
- => Battery charging system
- => Certified by KFDA, CE, FDA

Component hand piece hand piece is a user-friendly, Automated electronic with rechargeable battery, which allows it to lose the wire. It is a semi-permanent and both of needle tips and the guides are interchangeable.

#### Technical Specification

- => Dimensions hand piece-177.5mm x 34.88mm x 34.88mm  
cradle-R0.65mm x80.65mmx68.21mm
- => Weight : 530g (including package)
- => Battery lithium polymer : 500mAh
- => Duration recharging time : 3hrs
- => Continuous using time : 3hrs
- => Performance high speed mode : 700 cycle/min  
Low speed mode : 412 cycle/min

#### Product include Guide

- => hand piece
- => 3 recharging cradle
- => 1 power code
- => 6 guides
- => 5 needle tips
- => Are interchangeable to adjust the length of the needles.
- => it depends on different areas or purposes of treatment.
- => it also keeps the right angle of the needle insertion.





**KRUPA**

# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com



**Auto-MTS Therapy** is Derma Rejuvenation Therapy. Auto-MTS is Non-surgical skin Rejuvenation treatment which utilizes Micro Rolling Needle. Change the future of derma rejuvenation Derma Pen is a solution for the meso-roller.

- ✓ Portable Wireless
- ✓ Automatic micro needle
- ✓ Therapy care system
- ✓ Low pain
- ✓ Depth control available
- ✓ Competitive price
- ✓ Rechargeable

## Specifications

- ❑ Multi-needling system stimulating skin
- ❑ Depth adjustable (0.25mm~2.0mm) only with 1 needle
- ❑ Much decreased pain compared to meso roller
- ❑ Excellent improving effect
- ❑ Battery charging system
- ❑ Certified by KFDA, CE, FDA, 3patents

## Benefits

- ❖ Skin rejuvenation
- ❖ Enhancement of skin elasticity
- ❖ Wrinkle reduction
- ❖ Acne scar reduction
- ❖ Stretch mark reduction
- ❖ Stimulation of hair growth
- ❖ Whitening
- ❖ Administration of active ingredients

## Technical Specification

- ✓ Dimension hand piece-177.5mm x 34.88mm x 34.88mm cradle-80.65mm x 80.65mm x 68.21mm
- ✓ Weight : 530g (including package)
- ✓ Battery lithium polymer : 500mAh
- ✓ Duration recharging time : 3hrs
- ✓ Continuous using time : 3hrs
- ✓ Performance high speed mode : 700 cycle/min Low
- ✓ Speed mode : 412 cycle/min

## Parameters

- ◆ Adapter : 6V - 300mA
- ◆ Needle Size : 0.25~2.50mm adjust
- ◆ 2pcs Lithium battery
- ◆ Size : 14.7x 3.1cm
- ◆ Package : 18.5x15.5x5.5cm

## Specifications

1. Wrinkle removal, acne removal, improve scar, rejuvenation
2. For the narrow focus of the site can be treated(glabella, depth, wrinkle, scar)
3. Shrink pores
4. treatment of stria gravidarum
5. Lifting
6. Treatment of acne scars and healing wounds

## Technical Specifications

- Voltage: 5V=1000mA, 12VA
- Input: AC100-240V, 50-60Hz
- Size: 20\*20\*140mm
- Weight: 555g, pen: 82g
- Micro needles length: 0.25,0.5,0.75,1.0mm
- Charge Time: 40 Minutes
- Battery Life: 4 Hours
- Battery Capacity: 750mh
- Micro Needle: 12 Pin
- Accessories: Machine, Plug, Battery, Charger

## Guides

- ❖ Are interchangeable to adjust the length of the needles.
- ❖ It depends on different areas or purposes of treatment.
- ❖ It also keeps the right angle of the needle insertion.



ISO 13485 (Medical Standard)

Certified By :



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Derma Roller



Derma Roller is the world's leading Micro-Needling is non-invasive, Safe, Natural Cosmetic procedure. It activates stem cell to creat new/Pump Healthy skin through a process called Induction Therapy.



### Advantages

- ✓ Minimal risks, Safe, No risk of burns
- ✓ Cost effective, Less then laser treatments
- ✓ Short healing period
- ✓ Performed under topical anesthesia
- ✓ No permanent skin damage
- ✓ No increase in sun sensitivity
- ✓ Thin or pre-lasered skin can be treated
- ✓ Suitable for all type of skin
- ✓ All areas of the Face, Neck & Body can be treated



### 540 Needle

- ✓ Material of needle : Titanium
- ✓ Material of handle/Roller : PC+ABS
- ✓ Quantity of needle: 540pcs
- ✓ Approved by: CE & FDA
- ✓ High sealing sterilization packing

Damaged Skin  
Scars, Fine lines,  
Wrinkles, Stretch  
marks, Cellulite

Micro Needles  
Stimulate the skin  
to produce collagen


The new collagen  
produced starts moving  
upwards towards  
the skin surface

Collagen production has  
increase and is evenly spread  
in skin resulting to a smoother  
texture, small pore size,  
reduced scars, fine lines  
and wrinkles.

### 200 Needle

- ✓ Material of needle : Titanium
- ✓ Material of handle/Roller : PC+ABS
- ✓ Quantity of needle: 200pcs
- ✓ Approved by: CE & FDA
- ✓ High sealing sterilization packing

#### MOVING DIRECTION OF DERMA ROLLER



Before



After



#### APPLICATIONS IN CLINIC

	0.25mm	0.5mm	0.75mm	1.0mm	1.25mm	2.0mm	2.5mm
Cell Regeneration	✓	✓	✓	✓	✓	✓	✓
Anti Age / Anti Wrinkle	✓	✓	✓	✓	✓		
Stretch marks				✓	✓		
Hyperpigmentation					✓		
Acne Scar Light	✓	✓	✓	✓	✓		
Acne Scar Deep				✓	✓		
Hair loss	✓						
Hair Regrowth	✓						
Cellulite						✓	✓

AN ISO 9001:2008 CERTIFIED COMPANY

305.panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Derma Scope


✓ LED Light	: 8
✓ Wheel	: Focus wheel
✓ Switch	: Led light On/Off
✓ Button	: For Snapshot
✓ Image Sensor	: 2.0 Megapixels (Interpolated 2M)
✓ Still Capture Resolution	: 1600x1200, 1280x960, 800x640x480, 352x288, 320x240, 160x120
✓ Color	: 24bit RGB
✓ Flick Frequency	: 50/60 Hz
✓ Frame Rate	: Max 3 of/s under 600 Lus Brightness
✓ Magnification Ratio	: 40x to 800x
✓ Shutter speed	: 1Sec to 1/1000 sec
✓ Video Format	: AVI
✓ Photo Format	: Jpeg or BMP
✓ White Balance	: Auto
✓ Exposure	: Auto
✓ Light Source	: 8 LED's
✓ PC Interface	: USB 2.0
✓ Power Supply	: 5v Dc from USB port
✓ Supply Window	: 7/Vista/XP 112x30x30mm(approx)



Krupa's Derma Scope is a full fledged CCD Camera to view skin & Hair,Scalp condition etc. At a very high level of magnification & Resolution. The images through a vertical optical zoom lens which two heads gives range of magnification level 50x & 200x live zoom. It has a built in light source for illumination. The scope is extremely light weight & User friendly during the procedure.

This machine can be attached to computer/Laptop,LED Screen,Monitor keep record of pre-post photos & microscopic images.

### Derma Scop images



The most imaging High quality technology for skin analysis. It's the smart choice for your Remote Clinic.

**KRUPA MEDI SCAN**

AN ISO 9001:2008 CERTIFIED COMPANY

305,panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

# Derma Touch




**Features**

- High speed tower model with mobile unit
- Virtually no side effects
- Effective on all skin color & types
- No anesthesia or recovery time required.
- Immediate result & continuing improvement
- Patient return immediately to daily activities
- can be performed quickly (about 15 to 20 Minutes)
- 2 in 1 model- Crystal+Diamond
- 3 different outlets 100/140/180 microns
- Non blocking systems
- Automatic crystal flow

**Applications**

- Active acne scar revision
- Hyperpigmentation
- Stretch marks
- Dull, Congested, Thick or Blemished skins
- General Dermabrasion
- Face & body peeling lined
- Wrinkled or sundamaged skin
- Tattoo removal
- Age spot from the face
- Remove black & white heads
- Remove Un even skin tone
- Treat melasma
- Remove burns or post truma scars.

**Technical Features**

Hand Piece Material	As per medical grade.
Hand Piece	Pencil Type.
Tips Material	Standard Medical Grade ( can be sterilized ).
Diamond Probe	Pencil Type.
Diamond Tips	3 Pcs as per medical standard.
Dihumidifier	With Hydro Beads.
Delivery Angle	90 degree.
Exfoliation Agent	Pure Aluminum Oxide Sterilized Crystals.
Particle Distribution	100 / 140 / 180 Microns.
Containers	Fresh Supply ( 3 ) Jars. Waste Collection jar ( 1 ).
Vacuum Pump	2 nos imported Noise free pump.
Powder Flow	Automatic as per Vacuum Adjustment.
Vacuum	26 inch above.
Operating System	Micro Processor Base.
Electrical Supply	S.M.P.S 220V 50HZ. ( In Build Spike Guard Protection ).
Dimensions ( cm )	L 36 x W 30 x H 26.
Weight	18kg Approx.

Certified By :   ISO 13485 (Medical Standard)

## Digital Iontophoresis

Iontophoresis treatment are recognized across the world by dermatologist as one of the few successful hyperhydrosis treatment of Armpits,Hand & Feet,No Drugs,No Surgery.



Iontophoresis therapy for excessive sweating involves passing a low electric current through the affected body part(usually Hand & Feet).Iontophoresis is generally used for drug delivery.It is unclear why iontophoresis work for hyperhydrosis.One explanation is that minerals in the water get into skin and block the sweat ducts.

### Application

- ★ Hyperhydrosis using plain tap water
- ★ Drug delivery
- ★ Cutaneous Anaesthesia

### Accessories

- ★ Plastic for Iontophoresis (2Psc)
- ★ Electrode plates

### Technical Features

- ◆ Power - AC 220/240v 50Hz
- ◆ Microprocessor base operating system
- ◆ Digital procedure count down timer
- ◆ Intensity range: 0 to 30mA
- ◆ Spike guard protection
- ◆ Auto cut off facility with beep sound
- ◆ Digital force display
- ◆ Variable force adjustment

## Digital Radio Frequency 4.1MHz With Hair Removal



### Applications

- Moles,Warts,Skin tags
- Compound nevus
- Epidermal nevus
- Depressed scars
- Adenokeratosis
- Granuloma
- Birthmark
- Rynophyma
- Xanthelasma

Electrode set



### Accessories

- |                    |                               |
|--------------------|-------------------------------|
| • Probe            | : Imported pencil type for RF |
| • Electrocode set  | : Tweezer type for HR         |
| • Conductivity gel | : Straight needle             |
|                    | : for Hair Removal            |
|                    | : Round loop 5mm              |
|                    | : Round loop 8mm              |
|                    | : 45 degree needle            |
|                    | : Triangle loop 5mm           |
|                    | : Diamond loop 8mm            |
|                    | : Easy to operate             |
|                    | : For patient safety          |
|                    | : 3 pin Grounded plug         |
- Foot switch
  - Earthing plate
  - Mains code

### Technical Features

- Micro processor base operating system
- Pressure less cutting
- Minimum literal heat
- Micro smooth incision
- Reduced operating time
- Frequency: 4.1MHz
- Line voltage:220/240v AC/50Hz
- Power output: 130watt
- Mode of function:
  - 1.Coga: 60%coag-40%cut
  - 2.Cut : 70%cut-30%coag
  - 3.Fulguration : Available
- Probe: Imported Pencil type for RF  
Tweezer type for HR
- Foot switch activated
- Autoclavable hand piece
- Autoclavable Electodes
- Patient isolation transformer for complete safety
- No need of ground contact or skin contact of antenna plate
- Weight: 3.5kg

### Benefits for Hair Removal

- Permanent Hair Removal
- Time & Money saving
- Portable free from Waxing,Shaving,Lotions,Laser treatments,etc

### Description

The Krupa Hair Removal system is the most advanced hair removal system on the market. It removes unwanted hair for good, for the first time ever it is now possible to perform professional electrolysis professional results can be achieved using the Krupa hair removal system without the use of a needle. Tweezers are used instead to make the procedure safer and less painful clinical studies have shown that the Krupa hair removal system has equivalent effectiveness as professional needle electrolysis.

Before



After





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305.Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com



## Silectron Diode Laser



### Silent Features :

- 808nm \* The Brand Of "Dilas" Laser Bar Made by Germany.
- Fats-Safe-Reliable-Permanent Hair Removal
- cooling system ( Water + Air + Semiconductor + TEC) made By Italy
- Highly Durable Power Supply The brand Of "Meanwell" Taiwan
- Front Side American Imported "plug & Play" connector
- 8.4" Touch Screen Display With easy And User Friendly Operating Menu
- Big Spot Size. 12 X 20 MM cover treatment area instantly.

### Technical Specification

Model	: Ultra Cold
Input Power	: 2000W
Laser Type	: Diode Laser
Wavelength	: 808nm/810nm
Laser Output	: 800W
Pulse width	: 5-400ms
Energy density	: 1-120J/cm <sup>2</sup>
Spot size	: 12X12mm or 12X18mm
Repetition frequency:	: 1-20Hz
Cooling	: Sapphire touch cooling
Display	: 8.4 inch TFT chromatic screen
Cooling	: Water + Air + Semiconductor
	: Compressor Refrigeration cooling (Only in our Company)
	: TEC Cooling (Only in our Company)
	: Automotive Cooling blanket. (Only in our Company)





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

ISO,Photovoice certification, Air Management process, Safe Operation with respect to  
Australia road, Industrial law, Administration, ISO 9001:2008  
http://www.krupa.com / www.krupa.net / www.krupa.com  
E-mail: krupamediscan.com

## Diode Laser

## Suprimo

### Silent Features:

- World's First Diode Laser Hair Removal System.
- Fully-Auto-Heated-Pulsehead Hand Piece.
- Cooling system (Water + Air + Semiconductor) made by USA
- High Quality Power Supply The heart of "Mannat" Name
- Best Auto-Analyzer Imported "Plug & Play" connector
- 100% Original Design With easy And User Friendly Operating Menu
- Big Spot Size 12 X 20 MM gives Treatment area easily
- Reverse Three mode: Hair Removal / Skin Rejuvenation / Fast Hair Removal

Theory Treatment Principle:Diode laser removal technology is based on the selective dynamics of light and heat. The laser goes through the skin surface to reach the root of hair follicle; light can be absorbed and converted into heat damaged hair follicle roots; so that hair loss regeneration without injury surrounding tissue. It provides less pain, less operation, the safest technology for permanent hair removal now Diode laser machine adopt 808 nm, particular effectiveness hair melanocytes without injury surrounding tissue. The laser light can be absorbed by the melanin in hair follicles, and converted into heat, thus increasing the hair follicle temperature. When the temperature rises high enough to irreversibly damage the hair follicle structure, the hair will disappear after a period of natural physiological of hair follicles and thus achieve the purpose of permanent hair removal.



### Specification:

Laser Type	Diode Laser
Laser Wavelength	- 808nm
Out Put Power	- 2000watt
Pulse Width	- 5-500ms
Display Size	- 17" Touch Screen
Energy Density	- 1-10J/cm2 Adjustable
Spot Size	- 12 X 20 mm
Frequency	- 1-10KHz/Sec
Cooling	- Water + Air + Semiconductor
Temperature Of Probe	- 0-4°C
Weight	- 80Kg
Operating Voltage	- 230V / 50 Hz
Size	- (330X277X200)





## *Diode Laser*

Suprimo

#### **Editor's Feature**

- **Wetson®** The Brand Of "Highly™" Low Temperature In USA.
  - **PCB Anti-Pollution Protection With Thermostatic Control system (Water-air + Biomass) made By USA.**
  - **Highly™ Domestic Power Supply The Brand Of "Wetson®" Taiwan.**
  - **Black Box American Product "Plug & Play" connection.**
  - **PCB™ Smart Screen Display With Safety And User Friendly Operating Menu.**
  - **Big Input Size: 12V 20A With over protection area instantly.**
  - **Separate Three mode : Hot Thermostatic / Sun Regulation / Power Max Protection.**

Theory Treatment Principle: Sclare laser hair removal technology is based on the selective dynamics of light and heat. The laser goes through the skin surface. To reach the root of hair follicle, light can be absorbed and converted into heat damaged hair follicle tissue, so that hair loss regeneration without injury surrounding tissue. If presented this pain, easy operation, the safest, technology for permanent hair removal. sclare laser machines adopt 808nm, guarantee effective to hair melanocytes without injury surrounding tissue. The laser light can be absorbed by the melanin in hair follicles, and converted into heat, thus destroying the hair melanopigmentation. When the temperature rises high enough to irreversibly damage the hair follicle structure, the hair will disappear after a period of natural physiological of hair follicles and thus achieves the purpose of permanent hair removal.

## Specifications

Laser Type	Q-Switched Nd:YAG Laser
Laser Wavelength	532nm
Out Put Power	200mJ/pulse
Pulse Width	3-500ns
Display Screen	12" Touch Screen
Energy Density	1-100 J/cm <sup>2</sup> Adjustable
Spot Size	12 X 20 mm
Frequency	1-10Hz/Sec
Cooling	Water + Air + Semiconductor
Temperature Of Probe	0-4°C
Weight	60kg
Operating Voltage	220V / 50 Hz
Size	500E121120mm



Page 1

三



1





The advance Technology 808nm Diode laser. It is effective to absorb melanin and goes through dermis to reach to adipose tissue acting on different area and depth of hair. Excellent treatment results can be expected at the first treatment and suitable to all types of hair.

## Suprino Plus





### Smart Features

- 808nm - The Best 808nm Diode Laser Made by USA
- Auto-Duty Available Photorejuvenation Hair Removal
- Cooling system | Water + Air + Semiconductor made by USA
- Highly Durability Power Supply The brand Of "Mitsubishi" Taiwan
- Super Safe American Imported "Easy & Play" connector
- 10.4" Touch Screen Operation With Easy And User Friendly Operating Menu
- Big Area Size, 12x12 mm cover treatment area instantly
- Separate Three mode - Hair Removal / Skin Rejuvenation / Fast Hair Removal

Before

After



### Specifications:

Laser Type	Diode Laser
Laser Wavelength	808nm
Out Put Power	2000mW
Pulse Width	1-400ns
Display Size	10.4" Touch Screen
Energy Density	1-100 J/cm2 Adjustable
Beam Size	12 X12 mm
Photorejuvenation	1-10Hz
Cooling	Water + Air + Semiconductor
Temperature	50-55°C
Weight	10kg
Operating Voltage	100-120V



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

KRUPA

305,panorama complex, Nr Navjivan press, B/h Gujrat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8466555111 / 84666000006  
E-Mail: kms@krupamediscan.com

## Fractional Co<sub>2</sub> Laser

Ultra pulse treatment is a "Fractional treatment" that produce thousands of deep,tiny columns in your skin called microscopic treatment zones. Fractional in that it intensively treats old epidermal pigment cells without affecting any of the surrounding tissue,which reduces healing time about 24 hours and promotes the skin's natural process for creating new,healthier tissue to replace the damaged cells. Ultra pulse treatment is safe, non-procedure that allows you to rediscover the fresh healthy skin of your youth.

### Advantages

- Minimum space requirement
- Big 8.4" color LCD touch screen
- Humanized software control,Stable laser output
- Individual scan pattern for different applications
- Beam transport device with 7 articular arm
- Water+Air cooling system
- According to application select mode  
=>Normal mode =>Fractional mode
- Enlarge and reduce scanning graphics
- 6 group of parameter settings
- 3 types of scan methods  
1)free scan 2)Sequence scan 3)Max distance scan
- Exchange 3 kind of methods
- 10600nm wavelength protective laser glasses provided for the operator.

### Specifications

◆ Wavelength	: 10.600nm
◆ Controller	: 8.4 color touch LCD screen
◆ Pulsed Radio Frequency	: 0.530w
◆ Average state	: Ultra pulse mode
◆ Average power	: 30w
◆ Scan pattern	: Square,Rectangle,Round triangle,Oval,Diamond Shape,Line or Customized graphics
◆ Scan pattern size	: 0.1x1.1mm - 20x20mm
◆ Scan mode	: Free scan,Sequence scan,Maximum distance scan
◆ Spot size	: 0.2mm at the focus
◆ Pulse duration	: 0.1 - 10ms adjustable
◆ Spot distance	: 0.1 - 2mm adjustable
◆ Pulse energy (power)	: 1mJ to 100mJ adjustable on LCD screen
◆ Interval b/w scanning	: 0-6s
◆ Laser apparatus	: Sealed off laser device stimulated by direct current
◆ Condenser focus	: 50mm
◆ Beam divergence angle	: 0.3 mrad
◆ Aiming beam	: 635nm infrared ray
◆ Beam transport device	: 7 Articular arm
◆ Cooling system	: Water + Air
◆ Dimension of machine	: 53 x 57 x 138cm
◆ Weight	: 55kg
◆ Power input	: 220v/50Hz



### Application

- Acne scars
- Actinic keratosis
- Laser skin resurfacing
- Improve texture tone & pore size
- Smooth wrinkles around the eyes & mouth
- Rejuvenate your neck,chest & blemishes
- Reduction of age spot & Blemishes
- Remove intractable chloasmas & pigment
- Sun damage recovery



# KRUPA MEDI SCAN

An R&D perspective of the new economy

2006, разработаны стратегии для дальнейшего привлечения иностранных инвесторов в регион. Важнейшим из них является создание инновационного центра.

## Laser Hair Regrowth Machine

## Introduction of Laser hair Regrowth Machine

- 2000 plus different road losses; 90% from road traffic
  - Priorities: pedestrian crossings, children's areas
  - Separate cycle paths required especially along busy roads
  - Local road prioritisation - Children 100 (prioritise cycling)
  - Cycle right approach
  - No priority, no room either - just no consequences

第十一章

- Without pain treatment it diminishes your image and limits your confidence.
  - Increases blood flow which in turn increases hair growth in follicles not capable of producing hair.
  - Doctor has to tell 80% of patients requires additional pain activity for people with density hair in oily scalp.
  - Report suggests that poorly designed lasers increase incidence hair by up to 25% per hair producing much thicker and faster appearance.
  - Reduces oily scalp conditions.
  - Decreases recovery time after used in conjunction with hair transplants.
  - Safe, effective and permanent procedure.
  - No side reactions with under medical supervision.

#### Technical Specification

Color	670nm Red
Light Power	0-100mw Adjustable
Time	1-90 min Adjustable
Treatment panel	4 Pieces
Treatment Area	830mm * 200mm
Diode Laser Number	240Pcs super flux Japan
Gross Weight	40kg
Packing Size	100*49*51
Electrical	230v/50hz



INT'L PRACTITIONERS  
ASSOCIATION

CE SFDA



ISO 13485





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/H Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866551111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Diode laser Hair Growth device

### Clinical indication

With Slectron Diode Laser Treatment, there is increased by 73% of healthy hair and no side effect. Treatment Functions Our hair has its own life cycle and it only grows at its growth cycle. According to research, Slectron Diode Laser Photo-biological treatment is able to induce the hair follicle enter growth cycle. In addition, the hair that developed through stimulation is healthier and stronger than naturally grow hair. Clinic proves that only Slectron Diode Laser is capable of increasing the blood flow of scalp nutrient cell by 60%.

### Principle of the treatment

Laser light uses the scientific principle of photo-biotherapy. This process occurs when cells absorb light and repair themselves. Laser light penetrates deep into the tissues of the scalp, (it doesn't penetrate bone), stimulating microcirculation of blood supply, cell metabolism and protein syntheses. Clinical studies suggest that this cellular and sub-cellular activity breaks down the collection of DHT, helping to stop hair loss. Many long-term international clinical research programs have shown that bio-stimulation of 'cold beam' low-level laser gives numerous positive effects to the skin and hair. Research suggests that laser works by three mechanisms:

\*Laser light increases the energy available to cells so they take in nutrients faster and get rid of waste products . Not surprisingly, many scientists and physicians have been using low-level laser over the past 30 years to speed wound healing and produce new tissue.

\*Circulation, increasing blood and lymph circulation through vasodilatation. Microscopic studies have shown laser increases circulation and oxygenation of the blood to the scalp and hair bulb. Laser also removes calcification and blockages around the hair bulb. As well as increases cell replacement activity. Laser helps to improve hair in fullness , shine, body and elasticity. Therefore, problems such as excessively oily or dry scalp, dandruff and itchiness can also be reduced.

### Accessories

Skin Analyzer



Helmet



Medicine Sprayer



Microcurrent  
Massage Comb



### Technical Specification

Name	Diode Laser Hair Loss Therapy Device
Laser output power	650nm diode laser 320mw 808nm IR laser 2200mw
LED output power	630nm Red cool light 1000mw 460nm Blue Cool light 1000mw 580nm Yellow Cool light 2000mw
Treatment time	0-30min
Pulsing time range	1-10s
Pulsing Delay range	0-10s
Input power:	200w
LCD Screen	8inch multi color touch screen
Device Dimension	1630*550*500 mm
Package Dimension	1650*580*580 mm
Rated input voltage	AC230 10%,50Hz 1Hz
Net weight	35Kg
Package weight	55Kg

### Product Features

Lasers and PDT Lights

Complementary Medicine : Minoxidil

Easy To Operate

Comfortable Treatment Process

Remarkable Results

B  
E  
F  
O  
R  
E



A  
F  
T  
E  
R





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY  
305, Panorama complex, Nr Navjivan Press, B/H Gujarat vidhyapith  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Multipurpose Localized Phototherapy Unit Hand & Foot therapy



### Applications

- Psoriasis
- Vitiligo
- Mycosis Fungoides
- Atopic Dermatitis
- Lichen Planus
- Other skin disorders can be treated

### Features

- Less space requirement.
- Microprocessor base digital timer for accuracy
- UVA/NBUVB hour meter for all panels
- Individual On/Off switch for easy panels
- Imported mirror type reflectors ensure max irradiation
- Spike guard protection
- MCB ensure max safety for the patient
- MCB 10Amp
- Operating voltage: AC 230v/50Hz
- Power Consumption: 2 KVA
- No of lamps: 12 lamps
- UVA lamps: Philips
- The unit is mounted on a swivel castor wheels & locks are provided in the sides for easy operation
- Specially designed patient data cards are provided along with the unit
- Dimension-a) Height: 1500mm  
b) Width: 950mm  
c) Depth: 700mm  
d) Weight: 70kg Approx
- Protective goggles: 2Pcs
- High quality electrode chokes for max life of lamps
- Warranty 2 years of unit (not for Tubes).



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Handy NB-UVB Tube



### How to Use NB-UVB Tube

- ✓ Plug the A/C Socket in side of your tube.
- ✓ Remove the Upper Cover(Acrylic).
- ✓ Wear the protective goggles(Sun glass).
- ✓ Keep 1inch far from the part where you want to take the rays.
- ✓ Take the rays Max 10 min-Twice a day.
- ✓ Do not look direct to the NBUVB rays,it may damage your Eyes.
- ✓ No warranty of NB-UVB tube.
- ✓ Warranty only on Electronic parts for 1 Year.
- ✓ High Quality Electronics blast.
- ✓ Operating voltage: 230V/50Hz

### Indication:

- ★ Psoriasis
- ★ Vitiligo
- ★ Atopic Dermatitis

### Treatment Of Dose

Month	Time	Required Dose
1	5 Minute	Twice a day
2	6 Minute	Twice a day
3	7 Minute	Once a day
4	8 Minute	Once a day
5	9 Minute	One time alternate day
6	10 Minute	One time alternate day

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## IPL + RF System



12 X 40mm IPL PROBE



15 x 50mm IPL PROBE



TRI-POLAR RF PROBE



MONOPOLAR RF PROBE



### Hair Removal



Before



After

### Skin Rejuvenation



Before



After

### Pigmentation



Before



After

### Active Acne



Before



After

### Laser Treatment

- Remove Unwanted Hair
- Remove Pigmentation
- Non Surgical Face Lifting
- Remove Acne & Acne Spot
- Body shaping

- Skin Tightening
- Skin Rejuvenation
- Breast Enhancing
- Laser Photo Facial
- Remove Wrinkle

### Technical Specification

- 1 Light source : IPL
- 2 Fluence: 10-50 J/cm<sup>2</sup>
- 3 Pulse duration : 1-30MS
- 4 Pulse delay : 1-100ms
- 5 Pulse number: 1-15
- 6 Skin cooling : 0 - -4° C
- 7 Weight : 30 kg(net weight )
- 8 Output power : 2000w
- 9 Power : 230v/110v
- 10 Spot area: 12x40mm to 15x50mm  
+tripolarRF and monopolar RF
- 11 Output frequency-1-10 MHZ
- 12 Power – 1-100 J/cm<sup>2</sup>
- 13 Probe size – 16mm<sup>2</sup>, 32mm<sup>2</sup>, 36mm<sup>2</sup>, 40m<sup>2</sup>



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305.Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kcms@krupamediscan.com

## Laser Hair Comb



The Power Grow Comb is a breakthrough laser hair loss treatment that promotes the appearance of new hair growth. Using lasers and infrared LED lights, the Power Grow Comb helps stimulate hair follicle cells and improve scalp circulation.

In addition, vibration technology helps stimulate blood flow on the scalp and in turn produce fuller, healthier hair growth.

The comb itself has teeth that are positioned perfectly to create separation and allow laser LED to reach the scalp's surface and effect the hair follicles directly.

### Technical Data

- 7 Red 660 nanometer energy diodes
- 7 Blue 470 nanometer energy diodes
- High precision laser beam window
- Scalp stimulating vibrating action
- Set includes:

Power Grow Comb,  
Massage Comb,  
Stimulating Hair Brush,  
Manicure Set

### Direction

- Use the Power Grow Comb without applying gels or lotions to the hair first. Ideally use on dry hair.
- Use for 10-15 minutes every other day. Visible results take between 8-12 weeks to appear.
- Move the laser comb slowly around the head in the same direction of hair growth from front to back and top to bottom.



KRUPA

# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Localized UV Therapy Panel



### Application

- ❖ Psoriasis
- ❖ Vitiligo
- ❖ Lichenplanus
- ❖ Atopic dermatitis
- ❖ Micosisfungoides
- ❖ Other skin disorders can be treated

### Technical Data

- ❖ Small & portable dynamic model
- ❖ Microprocessor base digital timer
- ❖ Mirror type reflectors
- ❖ Power 220-240v 50Hz
- ❖ No of lamps 3
- ❖ Rotate up to 180 degree
- ❖ Easy to operate
- ❖ Auto cut off facility

## For Patient Use Only

### Application

- 1.Psoriasis
- 2.Vitiligo

- UVA~UVB Tube as per patient requirement
- 1Feet tube
- Operates on 220v/50Hz
- Protection cover attached
- long A/C main score

Certified By :



ISO 13485 (Medical Standard)



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Mesoderm

Enjoy Smooth,  
Soft and  
Glowing Skin



- Diamond microdermabrasion it is a treatment of mechanical peeling and moderate skin resurfacing. The treatment is performed using heads with diamond crystals of different levels of abrasiveness and stable adjustment of the vacuum. This is method of controlled removing dead cells of the epidermis layer by layer. Diamond microdermabrasion is used for problem skin resurfacing to remove acne, wrinkles, age spots, stretch marks, skin cleansing and prepare skin for further treatments by increasing skin absorption.
- This no-needle mesotherapy treatment including the transmission of active substances into the skin as a result of the phenomenon of electroporation. This mechanism uses TTS transdermal therapeutic system which consists in transmission particular substances into the skin by a layer of the epidermis and dermis using the optimal combination of sonotherapy (sonophoresis) with current pulses (iontophoresis).

### APPLICATION

- Acne, blackheads, clogged pores
- Hypersecretion of the sebaceous glands
- Enlarged pores
- Dry skin, desquamation
- Scars, striae, stretch marks
- Scars after acne
- Hyperpigmentation, depigmentation
- Skin aging, sagging
- Fine wrinkles elimination

### TECHNICAL FEATURES

- ✖ Digital version with nonblocking system
- ✖ Very portable and light weight
- ✖ SMPS power supply
- ✖ Inbuilt spike guard protection
- ✖ Microprocessor based procedure time setting
- ✖ Max vacuum 22 inch and above
- ✖ Power - AC 220/240v-50Hz
- ✖ Vacuum pump: Imported quality
- ✖ Function of mode: Dermabrasion & Meso Therapy
- ✖ Using angle: 90 degree
- ✖ Sterile Diamond tips: 3Pcs
- ✖ Ultrasonic Frequency: 1MHz
- ✖ Current range: 0.0mA - 30mA
- ✖ Meso Therapy Attachment: Face, Eye & Hair Probe
- ✖ Technology: Sonophoresis Transdermal Technology
- ✖ Weight: 3.5 kg approx



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

2005, Panchavati complex, Nr. Navrangpura junction, 380009 Gujarat, India.  
Ashram road, Navrangpura, Ahmedabad 380014.  
Ph: +91-9825000001 / +91-9825000002 / +91-9825000003  
E-mail: [kruvamediscan.com](mailto:kruvamediscan.com)

## Meso Therapy



No needle Mesotherapy are advance technology when using active ingredients as it will help infusion of the product deep into the dermal layer of the skin. Use Ultrasound vibrations at 300 million times per second. The Ultrasound waveform stimulates the skin tissue and improves the permeability of the cell membrane which also increases collagen production.

This treatment also improves the appearance of fine line & wrinkles. Ultra sound also generates frictional heat that will improve blood circulation. Helps to break down the fatty deposits to reduce cellulite.

### Features

- Best static ultrasound base Mesotherapy unit
- Most advanced technology for best result
- LCD display with feather touch keypad
- Dual frequency mode 1MHz & 3MHz
- Easy to operate
- No need for training
- 1MHz use for deep penetration
- 3MHz use for super facial & sensitive skin
- User friendly function of mode

### Technical Specification

Power input	: 200v/50Hz
Ultrasound frequency	: 1MHz to 3MHz
Display	: 30x4 LCD Backlit
Max output in cont mode	: 0w/cm <sup>2</sup> or 10w
Max output in pulse mode	: 0w/cm <sup>2</sup> or 21w
Unit key	: watts/cm <sup>2</sup> or watts
Pulse repetition frequency	: 10-180Hz in step of 10Hz
Pulse average ratio	: Duty cycle 0-100% in step of 5%
Treatment time	: 20 Minutes maximum
Transducer	: Ergonomic heads with 1 & 3MHz
Program	: Customer base designed protocols
End of therapy	: Buzzer

Certified By :



ISO 9001:2008 (Medical Standard)



KRUPA

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## + Dr. Meso Therapy Gun

Mesotherapy Meso gun originated from French medical doctor, Dr. Pistor is performed for obesity, plastic surgery and pain-care in many countries actively injecting small amount of medicine based on the theory that dermal layer, subcutaneous layer and muscle are same. Light weight and relatively small size make the hand easy to perform. Possible to use continuously since injecting device and control parts are separated. Practical and convenient cart makes performance easy at small room. Needle guider helps to control injecting depth. Possible to use any needle and syringe under 10cc. Micro process inside makes possible for various injection modes.

### Advantage

1. Light weight and relatively small size make the hand easy to perform.
2. Possible to use continuously since injecting device and control parts are separated.
3. Practical and convenient cart makes performance easy at small room.
4. Needle guider helps to control injecting depth.
5. Possible to use any needle and syringe under 10cc.
6. Micro process inside makes possible for various injection modes

### Main Features

1. The entire Process is completed without the need of surgery and anesthesia.
2. Does not cause the phenomenon of uneven skin
3. Does not cause the phenomenon of bleeding swelling and blood stasis
4. Non invasive treatment does not affect the normal work and life.

### Application

- \* Anti Wrinkle
- \* Skin Regeneration
- \* Striae Of pregnancy
- \* Bio Freckle Dispelling
- \* Bio Skin purify
- \* Bio Skin Whitening
- \* Body Sliming
- \* Skin Lifting.

### Item Specific

- Shape: Portable
- Catheter: 100pcs
- Power: 30W
- Packing size: 43 x 35 x 14cm
- Warranty: 1 year
- Voltage: 110V/220V
- G.W.: 5 kgs
- Product name: +Dr Mesotherapy Gun
- Color: White



### Technical Specification

* Power Output	30Watt
* Duty Cycle	For optimum Condition
* Mode	Superfacial layer
* Flow Rate	1-9 Grade
* Speed	1-6 Grade
* Mode	Point to point
* Flow Rate	10,20,40,60,80,100 ul
* Speed	1-5 Grade
* Input Power	230v/50hz



INTERNATIONALLY APPROVED

SFDA CE ISO 13485





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjeevan press, B/h Gujarat vidhyapeeth,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Micro Current



### Application

- ◆ Pigmentation
- ◆ Under Eye dark circles
- ◆ Oily skin
- ◆ Active acne
- ◆ skin lightening

### Technical Features

- Power - AC 220/240v 50Hz
- Microprocessor
- Power supply
- Digital procedure count down adjustable time
- Spike guard protection
- Digital force display
- Individual gel for application
- Digital force adjustment
- Digital change over switch

### Altered Pigmentation



Before

After

### Acne



Before

After

Certified By :



ISO 13485 (Medical Standard)



# KRUPA MEDI SCAN

KRUPA

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Ozone Therapy

The direct method of use procedure called sparkling which creates OZONE. ozone is a mild antiseptic & has a drying effect on the skin. This treatment would benefit the blemished skin & control acne.



### Use of Method

1. Indirect Method
2. Direct Method
3. Sparkle Method

### Technical Specification

- ✓ Micro Processor base Operating system
- ✓ SMPS Power supply
- ✓ Count down timer
- ✓ Adjustable force control
- ✓ Feather touch on/off switch
- ✓ Operating voltage 220-240volt 50Hz
- ✓ Frequency 10KHz
- ✓ Absorbed power-30vA
- ✓ Max output power - 40kv
- ✓ Output watt- 22w
- ✓ Operating temp- 10°- 50°
- ✓ Shock proof cabinet



### Treatments

1. Anti acne spot treatments
2. Wrinkle reduction
3. Anti aging treatments
4. Hair follicles stimulation
5. Oily skin treatment
6. Fade dark eye circles treatments
7. Dandruff treatments



# KRUPA MEDI SCAN

AN ISO 13485:2003 CERTIFIED COMPANY

305, Panchayat road, Virar East, Mumbai 401102, India  
+91-9820000777 / +91-9820000778  
E-mail: [kumar@krupamediscan.com](mailto:kumar@krupamediscan.com)

## Photo Dynamic Therapy



### RED Light

The light source with intense pulsed light (IPL) can be used in Photo dynamic therapy with red light. Red light has the ability to penetrate deep into the skin which stimulates melanin production & the production of prostaglandins that stimulates the inflammatory process. Leading to decrease in pigmentation.

### BLUE Light

The normal cellular oxygen deficiency which occurs are not too severe. If normal cellular oxygen deficiency which can be treated by Bluephoto & Blue light.

When the erythrocytes are excited by the intense blue light they become hyperactive enhanced effect is due to the P450 system. The heat is the heat of the electrons.

### Advantage

- Strong energy & density
- Bigger area irradiation are design of the light head. Suitable for the character of face.
- One equipments with two lights(red & blue) they key for two lights can be easily switched.
- No hurt to the skin cell.No side effect.
- No special nurse & normal dressing for skin after treatment.
- Fit for any kind of skin
- Simple operation no need training
- 100% safety
- Multi arrange of light source face treatment

### Treatment

- Removal : White, Melasma, Pigmentation
- Improvement : Coarse pores, Baby skin like tone, Thin wrinkles, Gray complexion etc.
- Dose : Fleck Sunburn, Age pigment etc
- Remove : Skin or pigment caused by High energy from SP, & laser treatment
- Remove : 1 & 2 grade acne

### Technical Parameters

- Light source : LED gene Biology light
- Output wave : RED-640nm = 450mw
- BLUE-470nm = 450mw
- Output power : 120mw/cm<sup>2</sup>
- Input voltage : 220v/50Hz
- Operate system : Microprocessor base
- Protection : Spike guard built



Photodynamic Treatment  
Before After



Certified By :  
ISO 13485 (Medical Standard)



## Portable 4S System

### E-Light + IPL + RF + Nd Yag Q-Switch Laser

#### Applications

- **E-Light**: Hair removal, Acne treatment, skin rejuvenation, Wrinkle removal, Pigmentation treatment, Vascular therapy, Breast liftup assistance, Laser photo facial
- **IPL**: big spot size for full body hair removal
- **RF**: Wrinkle removal, RF Photo facial, Skin rejuvenation, Anti aging treatment, Skin tightening, Body shaping, Weight loss
- **ND YAG LASER**: All colors Tattoo removal, Mole removal, Pigment, Melasma removal, Skin rejuvenation, Frackles removal, Acne spot removal, Birthmark removal, Carbon peel, Permanent Laser hair

#### Specification

##### ✓ **Basic Specification**

- Controller : 10.4" color touch LCD Screen
- Output Power : 2000w
- Language : Multi Software Language
- Cooling system : Semiconductor + Air + Water + Cooling Gel
- Skin cooling : <10° C
- Voltage : 220V / 110V

##### ✓ **E-Light (IPL+RF) System**

- IPL Energy : 1-50J
- RF Energy : 1-50J
- RF Frequency : 6MHz
- Pulse No : 1-6
- Handpiece number : 1pc with 5 Filters
- Wavelength : 480nm, 530nm, 590nm, 640nm, 690nm standard
- Spot size : 12x30mm<sup>2</sup>, 12x40mm<sup>2</sup> optional

##### ✓ **IPL System**

- IPL Energy : 1-50J
- Pulse No : 1-6
- Handpiece : 1pc
- Wavelength : 640nm Standard
- Spot Size : 15x50mm<sup>2</sup>
- Output Power : 1500w



E-Light  
12x40mm



Filters

##### ✓ **RF System**

- RF Energy : 0-50J
- RF Delay : 100-150ms
- RF Frequency : 6MHz
- Interval : 1-6s
- Handpiece No : 1pc with 2 probes
- Spot size : 16/32mm



##### ✓ **Nd.YagQ-Switch Laser**

- Laser Type : Nd.YagQ-Switch Laser
- Wavelength : 1064nm/532nm, 1064nm/10nm
- Indicator ofarming light : Unlocked ruby indicator
- Energy : 5000mJ



Spot Size      532nm      1064nm      10nm



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

3075, Juhu Seawoods Complex, Nr. Juhu Seawoods Phase 6th, Mumbai - 400089  
Address: Mumbai, Maharashtra, India - 400089  
Phone: +91 22 41426646 / 41426647 / 41426648  
E-Mail: KRUPAMEDISCAN@GMAIL.COM

## Q-Switch Nd-Yag Laser



### Treatment Range

Q-Switch laser can penetrate into deep layer of the skin so it is easily absorbed by Blue/Black & Green color pigment therefore it is very suitable to treat the pigment disease in the dermal layer for getting flat tattoos, Eyebrow tattoo, Moles & pigmented disease.

Q-Switch laser can penetrate into light layer of the skin therefore it is very easily absorbed by the Brown, Red & Deep color pigment.  
for Example: Eliminating freckles/Pigment spots & other light color blemishes.

### System Characteristics

- ✓ The Unique water cycle technology so it can be directly used without preheating.
- ✓ The micro computer control so it is very precise & safe.
- ✓ The plastic shell, so it has the feature of anti-static,anti-interference.
- ✓ The stable performance faster touch panel.
- ✓ The temperature & water protection, so it can give the protection of machine.
- ✓ The internal structure of modular design, so it is easily Tested/Repaired & Maintained.

### Technical Specification

Laser Type	Q-Switch Nd-Yag Laser
Laser Wave Length	1064nm/532nm
Pulse Energy	>800mJ / <800mJ
Pulse Width	10ns
Frequency	1.0Hz
Spot Diameter	1-6mm
Power Input	230v / 50Hz
Cooling System	Internal water + Air cooled
Net Weight	8kg
Energy	1000mJ

Certified By:



ISO 13485 (Medical Standard)



KRUPA MEDI SCAN

AN ISO 13485-2003 CERTIFIED COMPANY

CE MARKED MEDICAL EQUIPMENT MANUFACTURER & EXPORTER

# Q-switch Nd-Yag Laser



## Treatment Range

Nd-Yag laser can penetrate into the layer of the skin, so it is easily absorbed by black blood & deeper skin pigment therefore it is very suitable to treat the pigment disease in the skin layer for getting the better cosmetic effects like hyperpigmentation.

Nd-Yag laser can penetrate into the layer of the skin therefore it is very easily absorbed by the black skin & skin color diseases.

For example: Eliminating tattoo, Pigment spots & other skin color spots.

## System Characteristics

- The Unique water cooler technology to cool the cavity laser without pressurizing.
- The Micro computer control, so it is very precise & safe.
- The plastic shell so it has the features of Anti-static, Anti-infection.
- The stable performance higher output power.
- The temperature & water protection so it can give the protection of service.
- The internal structure of module design so it is easily cleaned/repaired & maintained.

## Specification

#	Laser type	Nd-Yag Q-Switch Laser
#	Wavelength	Double wavelength: 1064nm, 532nm
#	Indicator of coming light	Impacted ruby indicator
#	LCD controller	Key Press LCD display
#	Energy	1000mJ
#	Number of pulses	One
#	Frequency	1-10
#	Spot diameter	10mm
#	Cooling manner	Air + Water + Semiconductor
#	Program language	English
#	Dimension of machine	41x42x35cm
#	Dimension of package	41x42x35cm
#	Nett weight	10kg
#	Gross weight	17kg
#	Packing	Aluminum alloy case
#	Voltage	220V / 110V

## Accessories





## Q-Switch Nd-Yag Laser

Optional probe with Standard Accessories



1064 nm  
skin Rejuvenation

532nm

1064nm



Single Pulse 300mJ



single pulse 600mJ

### Treatment Range

1064nm laser can penetrate into deep layer of the skin, so it is easily absorbed by Blue, Black & Green color pigment therefore, it is very suitable to treat the pigment disease in the dermal layer for getting Rid tattoos, Eyebrow tattoos, Nevi & pigmented disease.

532nm laser can penetrate into light layer of the skin, therefore it is very easily absorbed by the Brown, Red & Deep color pigment.  
For Example: Eliminating freckles, Pigment spots & other light color tattoos.

### System Characteristics

- ▲ The Unique water cycle technology, so it can be directly used without preheating
- ▲ The Micro computer control, so it is very precise & safe
- ▲ The plastic shell, so it has the features of Anti-static, Anti-interference
- ▲ The stable performance father touch panel
- ▲ The temperature & water protection, so it can give the protection of machine

■ Laser type	: Nd Yag Q-Switch Laser
■ Wavelength	: double wavelength 1064nm, 532nm
■ Indicator of aiming light	: Imported ruby indicator
■ LCD controller	: 7.4" color touch LCD display
■ Energy	: 1600mJ
■ Width of pulses	: 6ns
■ Frequency	: 1-6
■ Spot diameter	: 1-8mm
■ Cooling manner	: Air + Water + Semiconductor
■ Program language	: English
■ Dimension of machine	: 37x50x36cm
■ Dimension of package	: 67x58x53cm
■ Net Weight	: 15kg
■ Gross Weight	: 29kg
■ Packing	: Aluminium alloy case
■ Voltage	: 220V / 110V

### Accessories





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

KRUPA

3005, Panoramic complex, Nr. Harghita street, 6th floor, Vasile Lupulescu, Andronache road, Andronache Tice, Arad, Romania - 3100014,  
Phone: +40244-066440 / +40244-055111 / +40244-055000  
E-mail: [kruvmedic@gmail.com](mailto:kruvmedic@gmail.com)

## Radio Frequency Cautery 2.8MHz



### Applications

- ✓ Moles, Warts, Skin tags
- ✓ Compound nevus
- ✓ Epidermal nevus
- ✓ Depressed scars
- ✓ Acne scarifications
- ✓ Granuloma
- ✓ Birthmark

### Technical Features

- ✓ Detachable switch for bi-polar
- ✓ Separate mode of Bi-polar & Monopolar
- ✓ Step switch intensity with light indicator
- ✓ Insulated electrodes
- ✓ portable & easy to operate
- ✓ Foot switch & hand switch activated
- ✓ Shock proof imported pencil type probe
- ✓ Foot switch activated with beep sound

### Features

- ✓ Pressure less-tetting
- ✓ Shock proof body
- ✓ Auto-detachable hand probe & electrodes
- ✓ Maintenance free technology
- ✓ No service required



### Technical Data

- ✓ Power AC: 220-240v/50Hz
- ✓ Frequency: 2.8MHz
- ✓ Operative with 100watt
- ✓ Function: Coag-Cut
- ✓ Model: Bipolar & Monopolar
- ✓ Intensity: step switch/light indication
- ✓ Hand switch: optional
- ✓ Patient isolation transformer for patient safety
- ✓ Size: 1.21cm X 1.6mm X 100cm
- ✓ Weight: 2.6 kg (Pencil)



Before



After



Before



After

Certified By:



ISO 13485 (Medical Standard)



### Digital Radio Frequency 3.5MHz With Hair Removal



## Technical Specification

Line input voltage	230v AC/50Hz
Frequency	5.5MHz
Power output	100watt
Electrode	Fully insulated
HIF probe	Imported pencil type
HIF probe	Imported twisted type
Couag	80% couag 40% oil
Cut	70% oil 30% couag
Polarization	Avalable

## Applications

- # Moles, Warts, Skin tags
  - # Compound Nevus
  - # Hypodermal nevus
  - # Depressed scars
  - # Achromic patches
  - # Granuloma
  - # Dermatofibrosis
  - # Syringomyelia
  - # Xanthelasma



## Technical Features

- 1) Non-invasive laser operating system
  - 2) Fully isolated Radio Frequency for patient safety
  - 3) Pressure less cutting for lateral skin damage
  - 4) No need of ground contact or skin contact on patient plane
  - 5) No need anaesthesia for normal Nose, Mouth or Lips
  - 6) Permanent hair removed by low temperature High frequency
  - 7) Effective on all types of skin & hair colors
  - 8) Painless Hair removal techniques
  - 9) Result oriented technology

## **Benefits for Hair Removal**

- ✓ Permanent Hair Removal
  - ✓ Time & Money saving
  - ✓ Portable free from Waxing, Shaving, Lotions, Laser treatments etc

• 100 •

The Auger Pen Network system is the most advanced and refined system of the optical automated armament for the gun. For the first time ever it is now possible to perform automated target performance results can be determined using the Auger Pen Network without the use of a media. Sensors are used instead to track the projectile's flight and test points. Previous studies have shown that the Auger Pen Network system is a reliable, efficient means for performance monitoring and analysis.





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY  
 305,panorama complex, Nr Navjeevan press, B/h Gujarat vidhyapith,  
 Ashram road, Income Tax, Ahmedabad-380014.  
 Ph: 9824406646 / 8866555111 / 8866000006  
 E-Mail: kms@krupamediscan.com

## Digital Radio Frequency 4.1MHz With Hair Removal



### Applications

- Moles,Warts,Skin tags
- Compound nevus
- Epidermal nevus
- Depressed scars
- Adenokeratosis
- Granuloma
- Birthmark
- Rynophyma
- Xanthelasma



### Accessories

- Probe : Imported pencil type for RF
- Electrode set : Tweezer type for HR
- Conductivity gel : Straight needle
- Foot switch : For Hair Removal
- Earthing plate : Round loop 5mm
- Mains code : Round loop 8mm
- 45 degree needle : 45 degree needle
- Triangle loop 5mm : Triangle loop 5mm
- Diamond loop 8mm : Diamond loop 8mm
- Easy to operate : Easy to operate
- For patient safety : For patient safety
- 3 pin Grounded plug : 3 pin Grounded plug

### Technical Features

- Micro processor base operating system
- Pressure less cutting
- Minimum literal heat
- Micro smooth incision
- Reduced operating time
- Frequency: 4.1MHz
- Line voltage:220/240v AC/50Hz
- Power output:130watt
- Mode of function:
  - 1.Coag: 60%coag-40%cut
  - 2.Cut : 70%cut-30%coag
  - 3.Fulguration : Available
- Probe: Imported Pencil type for RF  
Tweezer type for HR
- Foot switch activated
- Autoclavable hand piece
- Autoclavable Electrodes
- Patient isolation transformer for complete safety
- No need of ground contact or skin contact of antenna plate
- Weight: 3.5kg

### Benefits for Hair Removal

- Permanent Hair Removal
- Time & Money saving
- Portable free from Waxing,Shaving,Lotions,Laser treatments,etc

### Description

The Krupa Hair Removal system is the most advanced hair removal system on the market. It removes unwanted hair for good, for the first time ever it is now possible to perform professional electrolysis professional results can be achieved using the Krupa hair removal system without the use of a needle. Tweezers are used instead to make the procedure safer and less painful clinical studies have shown that the Krupa hair removal system has equivalent effectiveness as professional needle electrolysis.

Before



After





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

KRUPA

305, Panorama Complex, Nr. Navjivan Press, B/h Gujarat Vidhyapith,  
Ashram Road, Income Tax, Ahmedabad - 380 014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## New Latest Non Invasive Technology Tri-polar Rf & Multipolar Rf

No Downtime

Wrinkle Remove



Skin Rejuvenation



Cross fit



Laughing Line



Body Shaping



No Pain



No Surgery

Neck Eliminating



Breast Lifting



Skin Tightening



Inch Loss

### Functions

Monopolar Radio Frequency - 3MHz

Monopolar Radio Frequency Power - 50J

Bipolar Radio Frequency-3MHz

Bipolar Radio Frequency Power- 40w

Triplex Radio Frequency - 6MHz

Triplex Radio Frequency Power - 60w

Fourpolar Radio Frequency-6MHz

Fourpolar Radio Frequency Power - 60w

Multipolar Radio Frequency - 6MHz

Multipolar Radio Frequency Power -100w

### Application

Operating System : Micro Control Processor base

Power supply : S.M.P.S. 300Watt

Operating voltage : 230v/50Hz

Temperature Control : 45 Degree C Cutoff

Display : Digital

Timer : 1 to 99 Minutes

*Renew Your Confidence & Make your Life More Beautiful*

INTERNATIONALLY APPROVED

CE SFDA



ISO 13485





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

KRUPA

305, Panorama complex, Nr Navjivan press, B/h. Gujarat vidhyapith,  
Ashram Road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8886555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## RF-Excited Fractional Co<sub>2</sub> Laser

### Spiral Max

The New Developed "Spiral Max" Rf Excited Fractional co<sub>2</sub> Laser with(Metal Tube Made by USA).

Ultra pulse treatment is a fractional treatment that produces thousands of deep, tiny columns in your skin, called micro treatment zones. Ultra pulse is fractional in that it intensively treats old epidermal pigment cells without effecting any of the surrounding tissue, which reduces healing time (about 24 hours) and promotes the skin's natural process for creating newer and healthier tissue to replace the damaged tissue. Ultra pulse treatment is a safe, non-invasive procedure that allows you to rediscover the fresh, healthy skin for your beauty.

#### Application

- Wart, Viral Warts, Herpes Zoster, Nevus, Mole, Skin tags, Age spot
- Cutaneous Malignant tumour , Actinic Keratoses
- Eliminating acne scar and the scar by trasumulation or operation
- Remove Smooth wrinkle around the eyes and mouth
- Rejuvenate your neck, chest, back and hand etc.
- Laser Skin Resurfacing and rejuvenation
- Remove Intractable chloasma and pigment
- Improve skin texture tone and pore size.
- Remove underarm odour

#### Technical Specification

- Laser type : RF-Excited Fractional Co<sub>2</sub> laser (Metal Tube USA)
- Aiming beam : 635nm
- Beam delivery : 7-joint Articulated Arm
- Wave length : 10.600nm
- Display screen : 8.4 inch color Touch Screen
- Spot size : 0.12 - 1.2 mm
- laser power : 40Watt
- Scan size : 20 x 20 mm
- Scan pattern : Circle/Triangle/Square, Rectangle/Line/Ellipse.
- Cooling system : Air Cooling
- Operating mode : Ultra Pulse / Fractional
- Energy adjust : 1-100 mJ
- Pulse duration : 0.1-10ms
- Pulse repeat time: 1- 5000ms
- Overlap : 1- 20
- Pulse distance : 0.1 - 1.2mm
- Scan mode : Free Scan , Sequence Scan, Maximum distance Scan
- Power supply : 230volt / 50 hz
- Temperature : +5°C - +55°C
- Humidity : 30% - 80%
- Dimension : 55x82x125cm (LxWxH)
- Weight : 55Kg.



Scar Removal



Skin Rejuvenation

Result given by Dr. Nitin Patel, Mumbai, India.



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

KRUPA

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Rx-500



### Applications

- Super facial pigmentation
- Sun damaged skin
- Remove wrinkle on the face
- Remove acne scar
- Remove pigmented spot
- Skin & Body polishing
- Skin rejuvenation
- Remove enlarge pores



### Features

- Virtually no side effects
- Effective on all skin color & types
- No anesthesia or recovery time required.
- Immediate result & continuing improvement
- Patient return immediately to daily activities
- can be performed quickly (about 15 to 20 Minutes)

### Technical Features

- Digital version with nonblocking system
- Very portable and light weight
- SMPS power supply
- Inbuilt spike guard protection
- Microprocessor based procedure time setting
- Max vacuum 24 inch and above
- Power - AC 220/240v-50Hz
- Vacuum pump: Imported quality
- Crystal flow: Manual setting
- Function of mode: only crystal base
- Exfoliation agent: Sterile aluminum oxide crystal powder
- Delivery angle: 90 degree
- Sterile crystal probe tips: 3Pcs
- Dimension: H26 X W30 X L30 cm
- Weight: 7.5 kg approx

Certified By:



ISO 13485 (Medical Standard)



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Skin Analyzer



### Characteristics

- Imaging system: USA 1/3.25 inch CMOS
- Hardware: 2 mega pixel
- Maximum resolution: 1600\*1200
- Magnification factor: 50 times
- Operating temperature: 10-40
- Operating humidity: below 80%
- Power supply: USB 5V
- USB 2.0 port, plug and play, no need driving

Over the years, we have been engaged in offering a wide range of Skin Analyzer. The offered skin analyzer is distinguished in the market for its high strength and less maintenance. To accomplish extreme fulfillment, we sternly examine this analyzer on a chain of quality restrictions by our quality controllers.

Our skilled professionals provide these analyzers as per detailed specifications of our prestigious clients.

### Features

- => Support capture skin image under high purity visible light (IRGE 3) and polarized light
- => Analysis skin disorders in both epidermis layer and deeper layer
- => High definition image for operator to check skin issues including skin pore, pigmentation, vascular etc.
- => Automatically analyze skin issues as skin complexion, water ratio, oily ratio, wrinkles and texture, skin pores
- => Can simulate 3D (three-dimension) image of skin status
- => Customer information management, can show skin change trend with historical data
- => Integrate with cosmetic or treatment sales plan, can auto-generate prescription according to analysis result

### Specification

- => Image Captures System: High Quality CMOS Sensor
- => Image Definition: 1600 x 1200 Pixel
- => Visible Light Color: 7200 K
- => Polarized Light Support: Automatic analysis on
- => Analysis Item: Skin Complexion, Texture, Moisture, Oil, Pigmentation, Skin Pores
- => Interface: USB
- => Customer Management Software: Store Customer Info, Analysis Historical Data, Skin Trend
- => Cosmetic/Application Recommend: Auto Recommend Cosmetic/Application, with shopping cart
- => Weight: 135 g
- => Power Specification: 5 V/500 mA

Skin Analysis Report

Single Image Analysis

Professional Analysis

Capture Multiple Image

## Techno Beard Diode Laser

The Advanced Technology base 808nm Diode laser. It is effective to absorb melanin and goes through毛囊 to reach to adipose tissue, acting on different area and depth of hair. Excellent treatment result can be expected at the first treatment and suitable to all types of hair.



### Silent Features

- World's 1<sup>st</sup> Beard Diode Laser Pen Made by Germany
- Auto-Sterilization Function for Disinfection
- Auto-Main-Protection System (Water Flow + Water Level + Water Temperature)
- Cooling system (Water + Air + Semiconductor) made by Germany
- Unique Power Supply The Trend Of "Silent" Power.
- American Imported "Jaguar & May" controller
- 8.4" Touch Screen Display with easy And User Friendly Operating Menu
- Big Spot Size 12 X12 mm Cover Treatment area instantly.



### Specification

- |                        |                               |
|------------------------|-------------------------------|
| • Laser Type           | : Diode Laser                 |
| • Laser Wavelength     | : 808nm                       |
| • Output Power         | : 1000mW                      |
| • Display Size         | : 8.4" Touch Screen           |
| • Pulse Width          | : 10-1,400ms Adjustable       |
| • Energy               | : 1-120 J/cm2 Adjustable      |
| • Spot Size            | : 12X12 mm                    |
| • Frequency            | : 0.5-10Hz                    |
| • Cooling              | : Water + Air + Semiconductor |
| • Temperature Of Probe | : 0-3 C                       |
| • Dimension Of Machine | : 42*32*19cm                  |
| • Dimension Of Package | : 69*19*19cm                  |
| • Packaging Material   | : Aluminum Alloy Case         |
| • Weight               | : 30kg                        |
| • Operating Voltage    | : 220v / 50 Hz                |



# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

3005, Plot number 100, Sector 100, Panchkula, Haryana - 134101  
Phone: +91-172-2500000 / 2500001 / 2500002 / 2500003  
E-mail: [krupe@krupamediscan.com](mailto:krupe@krupamediscan.com)

## The Professional E-Advance E-Light System



### Features

- Imported Spare Parts & Korean Technology
- Five High Quality Filters for Different Treatments
- Contact cooling for patient comfort
- User friendly large touch screen
- Customize parameters settings possible
- Built-in preset parameters
- Portable & Result oriented
- High Quality Xenon lamp (USA)
- Large Spot area 12x40mm

### Technical Specification

- |                            |   |
|----------------------------|---|
| ➤ Light Source - E-light   | ➤ Display size - 7.4 inch (Colour Touch Screen) |
| ➤ Fluence - 1-500J/cm      | ➤ RF Energy - 1-50J                             |
| ➤ RF Frequency - 6MHz      | ➤ Dimension - 39x41x63cm                        |
| ➤ Pulse duration - 1-9.9ms | ➤ Weight - 25kg                                 |
| ➤ Pulse delay - 10-99ms    | ➤ Output power - 2000w                          |
| ➤ Pulse numbers - 1-8      | ➤ Power - 230v/50Hz                             |
| ➤ Skin cooling - 0-4c      | ➤ Lamp Quality - Xenon (USA)                    |
| ➤ Spot area - 12x40mm      |   |

### Filters

- 480 - 1200nm - Vascular Therapy
- 530 - 1200nm - Skin Rejuvenation
- 590 - 1200nm - Pigment Therapy
- 640 - 1200nm - Hair Removal
- 690 - 1200nm - Breast Lifting



INTERNATIONALLY  
APPROVED



SFDA

ISO 13485





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

**KRUPA**

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kcms@krupamediscan.com

## The Professional & Advance IPL System

### Advantage

- ❖ Five high quality filters for different treatments
- ❖ Contact cooling for patient comfort
- ❖ User friendly large touch screen
- ❖ Customize parameters settings possible
- ❖ Built-in preset parameters
- ❖ Portable & Result oriented
- ❖ High quality xenon lamp ( USA )
- ❖ Large spot area 12x40 mm
- ❖ Fast & Powerful cooling system
- ❖ LCD base total shots counter on the probe
- ❖ 2 trigger switch on the probe for easy operation
- ❖ Portable & Result oriented



### Technical Specification

- ❖ Light source - IPL
- ❖ Fluence - 1-50J/cm
- ❖ Pulse duration - 1-9.9 cm
- ❖ Pulse delay - 10 - 99 ms
- ❖ Pulse numbers - 1- 6
- ❖ Skin cooling - 0 - 4 c
- ❖ Spot area - 12x40 mm
- ❖ Display Size - 7.4 inch ( Colour Touch Screen)
- ❖ Dimension - 39x41x63 cm
- ❖ Weight - 25 kg
- ❖ Output power - 1500w
- ❖ Power - 230v/50Hz
- ❖ Lamp quality - Xenon(USA)

### Filters

- 480-1200nm - Vascular Therapy
- 530-1200nm - Skin Rejuvenation
- 590-1200nm - Pigment Therapy
- 640-1200nm - Hair Removal
- 690-1200nm - Breast Lifting



INTERNATIONALLY APPROVED



SFDA



ISO 13485





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr. Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Ultrasonic



### Application

- Face lift
- Breast lift
- Eye care
- Remove wrinkles
- Improve body lines
- Improves neck & cheek laxity

### Features

1. Spike guard inbuilt
2. SMPS power supply
3. Digital timer
4. Digital intensity
5. 3 types of pulse mode
6. Digital change over mode
7. Two probes (Face & Eye)
8. Powder coated shock proof body
9. easy to operate
10. Auto cut off facility
11. Feather touch operated

### Specification

- Input voltage
- Frequency
- Count down timer
- Intensity control
- Pulse mode

- : AC 220/50Hz
- : 3MHz output
- : 1-30 min
- : 0-9
- : I. Continuous
- : II. 1 Sec pulse
- : III. 2 Sec pulse



Before

After



ISO 13485 (Medical Standard)

Certified By:

Corporate office: Mr. Ravinder Singh, 10th Floor, Chancery Building, Sector 14,  
Ashoknagar, Gurugram, Haryana 1220014.  
Ph: 0124-4262444 / 22622211 / 22622200  
E-mail: [krupeamediscan@gmail.com](mailto:krupeamediscan@gmail.com)

## Whole Body Phototherapy Systems

Unit are available with  
UVB & NB-UVB  
combination therapy unit



### Features

- Low space requirement
- Separate switch for UVAB/NB-UVB panels
- Tubes are imported from Philips Holland
- UV cabinet made of high quality sheet with powder coated
- Microprocessor based digital timer for accuracy
- Independent timer master for each power
- Individual MCU switch for each panel
- Imported metal type reflectors ensure maximum irradiation
- Safe patient protection
- MCU (Main Circuit Breaker) ensures max safety for the patient
- MCU (Main Circuit Breaker) for patient safety
- UV protected glass can be provided for patient & operator
- Fully insulated control boxes
- Cooling fans provided
- Open side and front maximum ventilation
- High quality UV tubes for maximum life for tubes
- Remote control unit for easy operation
- Service & Repair by specially trained staff Engineers.

### Dimensions of unit:

	Height	Width	Depth
Outer Dimension	1900	1000	600
Inner Dimension	1900	900	500
Input power : 220V 110V			
Required power : 600W(16) 240 Volts 5A			
Weight : 100 Kg			

### Optional Features

- Whole body UV therapy systems are available with LCD Control Panel
- Computerized connectivity with software
- Whole body UV therapy are also available with square & Round shape
- UV Therapy with different combination (UVB & NB-UVB) can be provided as per your choice up to 48 tubes.

### Applications

- Psoriasis
- Vitiligo
- Mycosis Fungoides
- Atopic Dermatitis
- Parapsoriasis
- Lichen Planus





# KRUPA MEDI SCAN

AN ISO 9001:2008 CERTIFIED COMPANY

305, Panorama complex, Nr Navjivan press, B/h Gujarat vidhyapith,  
Ashram road, Income Tax, Ahmedabad-380014.  
Ph: 9824406646 / 8866555111 / 8866000006  
E-Mail: kms@krupamediscan.com

## Wood Lamp



### Diagnostic Uses

- ✓ Tinea versicolour
- ✓ Erythrasma
- ✓ Pigment alteration
- ✓ Vitiligo
- ✓ Leprosy

### Technical Specification

➢ Cabinet	: Molded(Acrylic)
➢ Flourescent tube	: 6 inch x 2Pcs
➢ Blue UV tube(360nm)	: 6 inch x 2Pcs
➢ Magnifying glass	: 65nm Diameter
➢ Choke	: Electronics
➢ Choke Capacity	: 8 watts x 2
➢ Tube Consumption	: 4 watts each
➢ Reflector	: Mirror type
➢ Size HxLxW	: 55nm x 160nm x 160nm
➢ Body	: Shock proof
➢ Weight	: 900 gram(Approx)
➢ Power input	: 170v-240v AC