

BPM-106 Digital Blood Pressure Monitor

Features:

- One button operation; just sit and press start for the blood pressure monitor to start.
- Very low noise during reading.
- ➤ User memory for 2 people; Storage for 2*120 readings.
- Arrhythmia Detection.
- Compatible with both batteries and adapters.
- Huge Screen to make it easy for everyone to read the blood pressure readings.
- Tempered glass to protect the screen.
- Ergononmic Design.
- Backlit Display.

How to Use:



Put the cuff on your arm cuff is in the correct direction.



Don't make it too tight or bend the tube



Make sure that the bottom end of the cuff is 2-3 cm from the elbow joint.



The pipe should be placed in the middle of the arm.



Place your arms flat on the table with your palms up and your body sitting straight. The center of the cuff should be placed at the same height as your heart position. Please note that the air tube on the cuff is not twisted.

Specifications:

Inflatable Device Measurement Mode Measurement Range

Accuracy

Resolution **Power Supply** Low Voltage Direction Storage Temperature and Moisture Using Temperature and Moisture

Automatic Inflation by internal pump Oscillometric method

Pressure: (30~290) mmHg; (4~37.3)kPa

Pulse: (40~199)beats/minute

Pressure: ±3mmHg

Pulse: ±5% of display reading

1mmHg

4* AA batteries

If there is the signal, replace the battery (-10~+55) °C Moisture: (10~95)% RH (5~+40) °C Moisture: (15~85)% RH

Contact us:



Mobile No./ WhatsApp No. : +91 9310113243





Kannu Impex (India) Pvt. Ltd.

