



One Button Operation



Heartbeat Waves



Oxygen Level



Easy to Use



Accurate Reading

## FTP-103 Fingertip Pulse Oximeter

### Features:

- One Button Operation; just sit and press start for the Fingertip pulse oximeter to start.
- Low Battery Indicator.
- Pulse Rate Waveform Display.
- Four Directions Mode.
- Large OLED Display.
- Low Power Consumption.

### Specifications:

Display Type	OLED Screen
Power Requirements	Two AAA alkaline batteries
SpO <sub>2</sub>	Measure range: 35-99% Resolution: 1% Accuracy: 70%- 100%. ±2% 35%-69% no definition
Pulse Rate	Measure range: 35-250 bpm Resolution: 1bpm Accuracy: ±2bpm Pulse intensity: bar-graph
Environment Requirements	Operation temperature: 5~ 40 °C Storage temperature: -10~ 50 °C Ambient humidity: 15-80%

### How to Use:



Open the pulse oximeter



Place your finger in the oximeter



Press the power button



Wait for the results.

### Contact us:



Mobile No./ WhatsApp No. : +91 9310113243



Web : [www.klifecare.com](http://www.klifecare.com)



Email : [info@klifecare.com](mailto:info@klifecare.com)

Kannu Impex (India) Pvt. Ltd.

412, Suncity Trade Tower, Sector-21, Gurugram-122015