

# NEB-102 Nebulizer

### Features:

- Strong built, easy to use and carry, all models with handles
- High performance nebulizers
- Fit for both adults and children
- Low noise production

## Specifications:

**Medication Capacity** Sound Level Nebulization Rate Flow Rate

**Operating Temperature** Operating Humidity

10 ~ 20 mL

≥ 60 dbA

0.2 mL/minute

3.5-6 Litre/minute

10°C ~ 40°C

10% ~ 95% RH

## How to Use:



the nebulizer Chamber.



Connect the tubing to machine and to the medicine chamber.



Close the cap of chamber attach the mask.



and use the nebulizer.



medication. Sit up and keep inhaling the medication until the mist stops.



wash the medicine cup and mouthpiece with warm soapy water and rinse them with water or towel.

### Contact us:



Mobile No./ WhatsApp No. : +91 9310113243



