

## NEB-103 Nebulizer

## Features:

- Strong built, easy to use and carry, all models with handles
- High performance nebulizers
- → Fit for both adults and children
- Low noise production

## How to Use:



Place the medication into the nebulizer Chamber.



Continue to inhale the

inhaling the medication

until the mist stops.

medication. Sit up and keep



Close the cap of chamber attach the mask.

Turn off the nebulizer and

wash the medicine cup and

water and rinse them with water or towel.

mouthpiece with warm soapy



Turn on the air compressor and use the nebulizer.

Contact us:

Mobile No./ WhatsApp No. : +91 9310113243

Web : www.klifecare.com

Email : info@klifecare.com

## Kannu Impex (India) Pvt. Ltd.

9 412, Suncity Trade Tower, Sector-21, Gurugram-122015

Specifications:

Medication Capacity Sound Level Nebulization Rate Flow Rate Operating Temperature Operating Humidity 10 ~ 20 mL ≥ 60 dbA 0.2 mL/minute 3.5-6 Litre/minute 10°C ~ 40°C 10% ~ 95% RH