









- 
Massage
- 
Cupping
- 
Acupuncture
- 
Tapping
- 
Scraping
- 
Available for
2 Users

TTD-101 TENS Therapy Device

Features:

- Five treatment modes
- Tailored pain therapy with 6 pain modes.
- Dual-channel treatment available for two people at the same time.
- Adjustable intensity of treatment at any time.
- Fall-off alert for electrode pad, followed by automatic suspension of treatment.
- Automatic shutdown in case of no-operation within 2 minutes except working state, more energy efficiency.
- Exquisite and compact appearance, light weight and portability
- Low battery indication.

Specifications:

Operating Mode	Continuous Operation
Treatment Time	5-90 Minutes
Intensity	0-15, Adjustable Interval; 1
Operating Voltage	DC 6V
Pulse Width	0.15-0.20ms, Accuracy ± 30%
Pulse Frequency	1-120Hz
Dimensions	121.8x6.6x22.2 (mm)
Electrode Pad	50x50 (mm)
Host Unit Weight	86g (Without Batteries)

How to Use:



Place the batteries



Press the start button.



Select the treatment area, and mode.

(Electrode Placement Chart)



Contact us:

-  Mobile No./ WhatsApp No. : +91 9310113243
-  Web : www.klifecare.com
-  Email : info@klifecare.com