

WONDFO BIOTECH
WeAreWorkingForYourHealth

www.rajbiosis.com
+91-998-3333-469

Wondfo

Guangzhou Wondfo Biotech Co., Ltd.
No. 8 Lizhishan Road, Science City,
Luogang District, 510663 Guangzhou, P.R.China
Tel: (+86) 400-830-8768
Website: en.wondfo.com.cn
E-mail: sales@wondfo.com.cn



Rev. A1 Rel.: 2022/08/15

Wondfo reserves all rights. All contents of this brochure are copyrighted by Wondfo.

Wondfo

Thyroid Function

Finecare™

Thyroid Function Rapid Quantitative Test

Regular screening of thyroid function!



Thyroid Function Rapid Quantitative Test

- ✓ Direct detection of Thyroid hormone levels
- ✓ Help diagnosis Graves' Disease, Hashimoto's Disease, Congenital hypothyroidism (CH)
- ✓ Support whole blood samples; Results ready in 15 minutes
- ✓ Room temperature storage up to 24 months

CLINICAL APPLICATION

- Diagnose and monitor thyroid diseases
- Find the cause of an underactive thyroid gland
- Keep track of thyroid function in hyperthyroidism treatment
- Double-check the diagnosis of congenital hypothyroidism


APPLICATION SCENARIOS



ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Reaction Time	Detection Range	Storage
T3	W231	Whole blood/ Serum/Plasma	75 µL	15 minutes	0.61~9.22 nmol/L (0.4~6.0 ng/mL)	4~30°C for 24 months
T4	W232				12.87~300 nmol/L	
TSH	W220				0.1~100 mIU/L	
fT3	W200				0.40-50.00 pmol/L or 0.26-32.55 pg/mL	
fT4	W259				1.00-100 pmol/L or 0.078-7.77 ng/dL	

TEST PROCEDURE



ADAPTATION INSTRUMENTS

1. Finecare™ FIA Meter Plus (FS-113)
2. Finecare™ II Plus SE(FS-114)
3. Finecare™ III Plus (FS-205)

