

BiPAP system specifications guide

System One sleep therapy platform

Philips Respironics revolutionized sleep therapy by introducing bi-level positive airway pressure technology to treat obstructive sleep apnea. Emulating natural breathing with pressure relief at inhalation and exhalation, this proven technology does more than provide highly effective therapy; it delivers a level of comfort that can mean the difference between compliance and noncompliance for patients having difficulty adapting to traditional CPAP therapy.

As part of the System One sleep therapy platform, the bi-level technology in our BiPAP Auto and BiPAP Pro devices adds an additional level of intelligence to the industry's most advanced system. And the System One

advancements continue with remote switching between bi-level and CPAP modes, increased humidity output and control, and dramatically reduced noise levels.



System One	Modes	Comfort features					
		Flex comfort	Flex demo modes	System One humidity control	System One resistance control	Ramp (5-45 minutes)	Altitude adjustment
BiPAP Pro	CPAP Fixed bi-level	Bi-Flex	BiPAP Bi-Flex	Standard	Standard	Adjustable	Automatic
BiPAP Auto	CPAP Fixed bi-level Auto bi-level	Bi-Flex	BiPAP Bi-Flex	Standard	Standard	Adjustable	Automatic

System One	Encore patient management					
	Data transfer modes		Data capture			
	Phone-in EncorePro 2	SD card EncorePro 2	Patterns of use	Advanced event detection and reporting	Advanced event response	Flow waveforms
BiPAP Pro	Hours of use	Standard	Standard	AHI, FL, RERA, snore, leak, CA, OA, HY, PB	N/A	Standard
BiPAP Auto	Hours of use	Standard	Standard	AHI, FL, RERA, snore, leak, CA, OA, HY, PB	FL, RERA, snore, leak, OA, HY	Standard

Legend

AHI: Apnea Hypopnea Index
CA: Clear Airway Apnea
FL: Flow Limitation
HY: Hypopnea
OA: Obstructed Airway Apnea

PB: Periodic Breathing
RERA: Respiratory Effort Related Arousal

PHILIPS

RESPIRONICS

sense and simplicity

Specifications		BiPAP Pro IN651S	BiPAP Auto IN751S
Pressure range:	4 to 25 cm H ₂ O	✓	✓
Flex pressure relief:	Bi-Flex	✓	✓
Ramp time:	0 to 45 min. (5-min. increments)	✓	✓
Starting ramp pressure:	4 to EPAP; patient-adjustable 4 to min. EPAP pressure	✓ N/A	✓ ✓
Dimensions:	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H	✓	✓
Weight:	3 lbs; 1.36 kg (without humidifier)	✓	✓
Filters:	Foam and optional ultra-fine	✓	✓
Device set-up:	LCD/control wheel/push button	✓	✓
Data storage capacity (minimum):	Display: 7- and 30-day averages SD Card: >1 year On-board: up to 6 months	✓	✓
Compliance meter:	Breathing detection	✓	✓
Altitude compensation:	Automatic	✓	✓
Electrical requirements:	100 – 240 VAC, 50/60 Hz	✓	✓
Humidification:	Integrates with System One heated humidifier and pass-over humidifier	✓	✓
Warranty:	Two years; replace with new	✓	✓
Miscellaneous:	Lighted LEDs	✓	✓

Ordering information			Ordering information	
Devices	Part number International	Australia	Accessories	Part number
BiPAP Pro with humidifier, SD card, Bi-Flex	IN651HS	AU651HS	Foam filter kit	1063091
BiPAP Pro with SD card, Bi-Flex	IN651S	AU651S	Ultra-fine filter kit	1063096
BiPAP Auto with humidifier, SD card, Bi-Flex	IN751HS	AU751HS	System One heated humidifier	1056209
BiPAP Auto with SD card, Bi-Flex	IN751S	AU751S	Data management accessories	Part number
Accessories	Part number		SD card (10 pack)	1063859
System One Performance Tube 15 (15 mm)	1074400		SD card mailer (10 pack)	1065146
System One Performance Tube 22 (22 mm)	1032907		SD card cover (attached to device)	1063858
Device side panel	1063784		SD card reader	1073915
Water chamber	1063785		Power cords (5ft, 1.5m)	Part number
Power supply	1063786		Power cord: Australia	1039012
Carrying case	1063857		Power cord: Europe	1039014
Power cord	1005894		Power cord: India	1039013
			Power cord: Japan	1039027
			Power cord: UK	1039026

Respironics, Bi-Flex, BiPAP, Encore, EncoreAnywhere and System One are trademarks of Respironics, Inc. and its affiliates. All rights reserved.



© 2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Broudy NA 11/10/10 MCI 4103802 PN 1068513

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9666 4444
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics