

#### **MANUFACTURES & EXPORTERS OF**

**ALL TYPES OF LAB & PHARMA EQUIPMENT** 

LAMINAR AIR FLOW, BIO SAFETY, FUME HOOD, DESPENSING BOOTH, PASS BOX, AIR SHOWER, WAX BLOCK, SLIDE CABINET & B.O.D. ETC.





# ORBITAL BOD SHAKER

### Orbital Shaking Incubator BOD Type





#### APPLICATIONS OF ORBITAL SHAKING INCUBATOR BOD TYPE

- •"ESCOLABS" brings the ultra-high-tech Triple walled Microprocessor controlled, solid-state Orbital Shaking Incubator down to earth.
- •Specially designed for Orbital Shaking under controlled Incubation / Temperature required for Tissue Cultures Fermentation Studies, & Ageing Tests., Life science applications & growth studies.

#### SALIENT FEATURES & CONSTRUCTION OF ORBITAL SHAKING INCUBATOR BOD TYPE

- •Robust construction outer chamber is made of Mild Steel duly pre treated finished with powder coated paint for lasting finish.
- •Inner chamber made of STAINLESS STEEL SHEET (304 Grade).
- •Fitted with durable coaxial blower for forced air circulation at triple walled back to maintain temperature uniformity in the chamber.
- •Outer front double walled door is insulated and is fitted with magnetic tape for air tight closing for no temperature loss, provided with lock and key arrangement.
- •Fitted with Universal Platform to accommodate interchangeable clamps of assorted sizes for different capacity of flasks..
- •Stationery shelves made of Stainless Steel (SS-304) for use of the unit as a standard B.O.D. Incubator, can be supplied at extra cost..
- •Door operated illumination lamp is fitted inside the chamber for easy visibility.
- •A full view inner plexi-glass door enables inspection and monitoring of inner chamber specimens without disturbing the process temperature.

۰H





High grade PUF/Glass Wool INSULATION between outer and inner chamber for minimal thermal losses.

- •Heating elements are placed in the path of moving air duly insulated from the body.
- •Control panel having all controls and circuitry are housed at the top fitted with Temperature controller cum Indicator Microprocessor based with LED display for Set Value (SV) & Process Value (PV) to select any desired temperature provided with Volt meter on the panel to read the incoming voltage.
- •Caster wheel mounted for easy portability.
- •Supplied with cord and plug.
- •Operating Voltage: 220 Volts AC (50 Hz).

# TECHNICAL SPECIFIC ACTIVES & ORDERING INFROMATION OF ORBITA SHAKING INCUBATOR BOTYPE

Usable Space in MM	600 X 600 X 600
Capacity	8 cu. Ft
Approx Volume (Ltr)	230
No. of Shelves	2
Temperature Range	5°C above ambient to 600 (Resolution 0°.C)
Temp.Controller	By Microprocessor Based PID Digital Temperature Indicator cum Contro
Shaking amplitude	25 mm.
RPM Display	Digital
Shaking Speed	Variable speed fro 28 RPM to 28 RPM
Temp. Accuracy	+1°C
Temp. Display:	LED Display for Set Value(SV) armodcess Value (PV)
Relay	Solid state electronic relay with protective heat sink.
Air Circulation	By forced convection system
Insulation	By High density PUF /Glass Wool
Operations	Nearly silent operation with ulltown vibration
Electric Supply	220/230V AC, 50/60Hz



Mob.: +91-9416985382, +91-7027572058, Email Id : sunilbrothers.sb@gmail.com, Website : www.sunilbrothers.co.in,





# **MANUFACTURES & EXPORTERS OF**



# SUNIL BROTHERS

(An ISO 9001:2015 Certified Company)

#4, SHANTI NAGAR, KARDHAN, AMBALA CANTT.- 133 001

Mob.: +91-9416985382, +91-7027572058

Email Id: sunilbrothers.sb@gmail.com, info@sbcompletesolution.com,

website: www.sbcompletelabsolution.com