



An




ISO 9001:2015  
CERTIFIED CO.

# MELTING POINT APPARATUS

The Melting Point Apparatus is a compact and efficient instrument used for determining the melting point of solid substances. It is widely used in laboratories for identifying purity and for pharmaceutical and chemical testing applications.



## ★ KEY FEATURES

-  Compact, lightweight and easy to operate
-  Uniform heating for accurate results
-  Clear visibility for precise observation
-  Suitable for capillary tubes
-  Durable construction for long service life
-  Ideal for laboratory and quality control applications

## 🏢 APPLICATIONS

  
Pharmaceutical  
Industries

  
Chemical  
Laboratories

  
Research &  
Development

  
Quality Control  
Testing

  
Educational  
Institutions

## ⚙️ SPECIFICATIONS


Application	Determination of melting point of solid substances
Temperature Range	Ambient to 300°C (Approx.)
Heating Medium	Silicone Oil / Glycerin / Paraffin Oil
Capacity	1 Test at a time
Capillary Tube	1.0 mm O.D. (Standard)
Thermometer	0°C to 360°C
Heater	Electric
Power Supply	230V AC, 50Hz
Construction	Powder Coated Body with Stainless Steel Top
Compliance	IS: 1448 (Part 10), IP Standards


## 📞 CONTACT US

-  [Info@eselindia.com](mailto:Info@eselindia.com)
-  [Sharma\\_electrical@yahoo.com](mailto:Sharma_electrical@yahoo.com)
-  [asinternational024@gmail.com](mailto:asinternational024@gmail.com)

 +91 9416183370  
+91 9416493069  
+91 9306996454

 International  
Business  
+91 7015448621

 62 Azad Nagar,  
Near Rampur Chowk;  
Behind Khalsa Tent House,  
Ambala Cantt,  
Haryana, 133001 India

 [www.eselinternational.org](http://www.eselinternational.org)  
[www.esellab.com](http://www.esellab.com)

also on [trade.gov.in](http://trade.gov.in)