



AICIL ROCKER SHAKER

A compact and light weight gel rocker finds extensive use in various laboratory testing procedures for a see-saw movement, and a angle tilt for gel staining and destining applications. suitable for bathing, blotting membranes, culture aeration etc. as well as molecular and other biological mixing applications. body made from powder coated CRC sheet duly powder coated with sea-saw platform furnished with a anti-skid platform rubber mat. with angle tilt of 4-5° and speed of 30/36 rpm.

Model No.	Platform Size
ARS-1	330x330 mm