



AICIL RH VIEW BOX

The RH View box , meets the RH slide test standards set by the American Association of Blood Banks. RH slide attains prescribed temperature, from 37°C to 39° C, within two minutes.

The view- area calibration range is from 45°C to 50°C. The Temperature is adjusted by a knob and measured on digital temperature meter accurately. The soft, glare free light provides optimal intensity for blood typing and cross matching. Light intensity is independent of the temperature control. The M.S Powder Coated Box is pivotally supportable in frame for gently rocking.

SPECIFICATION:

- Digital temperature indicator for easy accurate monitoring of viewing area.
- Viewing area temperature adjust easily to compensate for ambient temperature change (Range 45°C to 50°C).
- RH slide reaches temperature of 37°C to 39°C in two minutes, eliminating wasted heat up time.
- Soft, fluorescent, glare free light provides uniform illumination of slide content.
- View Box can be gently rocked using knob on support cradle.
- Sturdy, easy to-clean metallic powder compact case. Utilizes minimal bench space.

Model No.	ARHB-01
Light	Flourescent glare-free
Temp. Display	PID Digital temperature controller
Dimensions (LxWxH)	360x125x125 mm
Power	220V/50 Hz

FEATURES:

- Blood Typing.
- RH Determinations.
- Warming slides for Gram's stains.
- Tissue Typing.
- Reading Enzymes antibody screen in micro paltes.