



SEDIMENTATION SHAKER

AICIL SEDIMENTATION SHAKER

Specially designed for determination of sedimentation value of flour samples in order to measure protein quantity and quality based on swelling of protein globules as well as agglomeration of present starch into protein globules in a lactic acid medium. Fitted with digital RPM meter and digital timer. A cylinder carrier / Rack to hold 5 samples of glass cylinders are driven through a heavy motor. Price is without cylinder

Model No.	ASS-5
Main Body	CRC-PC
No. of Samples	To hold 5 Samples
Oscillation / min	40 oscillation /min
Tilt	30 degrees, each side
RPM indication	Digital RPM
Digital ON/OFF timer	To set upto 999 minutes
Supply	220/230 volts AC supply