## VARIANTS

Category	Variant name	Product code	Sponge size (mm)	Unit of sale
Cubes	Dental SL	GL111	10 x 10 x 10	32 Pcs/Pack
	ENT Size/ Dental XL	GL202007	20 x 20 x 07	10 Pcs/Pack
	Dental M	GL141407	14 x 14 x 07	10 Pcs/Pack
	Dental L	GL151515	15 x 15 x 15	10 Pcs/Pack
	Dialys	GL303010	30 x 30 x 10	10 Pcs/Pack
Regular	Standard	GL851	80 x 50 x 10	2 Pcs/Pack
	Export	GL751	70 x 50 x 10	2 Pcs/Pack
Custom Size	HEMOSPONGE is available in many variants other than mentioned above as per need of physician. Check our Master Catalogue or contact company consumer care for more details. Visit us at: www.goodwill-lifesciences.com			

# **EASY TO USE**



Clean the cavity



Place the sponge lightly with minimal pressure



Wait till haemostasis observed

# **ABOUT COMPANY**

We are manufacturer of Absorbable Gelatin Sponge USP with HEMOSPONGE brand. The company started its operation in 2016 at Surat, Gujarat, India. Manufacturing facility have well designed cleanrooms for production of sterile devices and also equipped with quality control/microbiological labs for quality check. Beside HEMOSPONGE dental cubes, we manufacture wide range of gelatin sponge to suit any surgical procedures. Visit our website for more details.

#### **GOODWILL LIFESCIENCES**

#### **Manufacturing Unit Address**

157-158, RJD Integrated Textile Park, Paal - Hazira Road, Ichhapore, Surat - 394 510, Gujarat, INDIA. Mfg Lic No.: G/25/2165

#### **Contact Details**

W: www.goodwill-lifesciences.com E: info@goodwill-lifesciences.com, hemosponge@gmail.com

O:+91-261-2974173, M:+91-9879744173





# HEMOSPONGE Dental Cubes

Resorbable Haemostatic Packing Sponge for Use in Dental Surgeries





# Goodwill Patient Safety First

# **FEATURES**



## Fast and Effective Haemostatic

Stops oozing blood in 2-4 minutes

## **Excellent Packing Material**

Provides physical supports and absorbs oozing blood from bleeding sockets

## **Highly Absorbent**

Absorbs 40 to 45 times the fluid/blood (1 gram of sponge absorbs 40 - 45 grams of blood).

#### **3D Matrix Structure**

UNIFORM POROSITY - better mechanical strength, superior absorbing power and better clot holding capacity

# Pure, Safe and Biocompatible

Made up of 100% Biocompatible material - extra pure parenteral grade **GELATIN ONLY** (Tissue collagen protein). complies biocompatibility standards - ISO 10993

# Non-toxic, Non-allergic and Inert Product

No tissue interaction,
Does not interfere natural tissue healing

# Resorbable (Absorbable)

Liquifies in 12-24 hours and drained slowly with salivation, 100% resorbable

# Biodegradable

When implanted deeper in tissue - metabolize by proteolytic enzymes within 3-4 weeks.

# Non-toxic Metabolites

The metabolites of degraded products are non-toxic and biocompatible

# **Keep Your Patient Smiling**

No need to remove after surgery - Provides higher patient comfort - Prevents secondary bleeding

## Neutral pH

Soak freely with any chemical, pharmaceutical or biologics for advanced effect

#### Use with Thrombin

Soak the sponge with thrombin or tranexamic acid for quick haemostasis. Soak with antibiotic – for local antiseptic effect

# Size and Shape

Wide range of size and shape to suit dental, ENT and other oral surgeries. (Available in 38 more sizes)

### **Dental SL**

Specially developed to fit bleeding sockets

## **Cost Effective**

10-15 INR per patient maximum

## Gamma Sterile, Ready to Use

No need to sterilize the product before use.

# **ACCELERATES NATURAL CLOTTING PROCESS**

Goodwill

Patient Safety First

**HEMOSPONGE** 

Pesorbable Haemostatic Packing Sponge for Use in Der

HEMOSPONO

HEMOSPONGE

HEMOSPONGE

Provides external surface to platelets to activate clotting process - no interaction with blood proteins