

VDRL SHAKER (LI-VDRLT-33 & LI-VDRL-33)

APPLICATION :

VDRL shaker is an essential equipment for a Venereal Disease Research Laboratory. It is used for the mixing of syphilis testing card, emulsion reaction, and complement fixation reaction while can also be used to blends of clinical diagnostic tests requiring gentle agitation.



FEATURES :

- Body made of MS duly powder coating.
- Shaker speed controller and timer (in LI-VDRLT-33 model only) provided.
- Carrier size 13" x 13".
- Fitted with geared motor with concentric motion.
- Spring clamp is fitted with the carrier.

TECHNICAL SPECIFICATIONS :

Model : LI-VDRL-33 & LI-VDRL-33	
Carrier Plate Size	13" x 13"
Shaker RPM	70-200 RPM
Timer	0-15 minutes (Timer facility on LI-VDRLT-33 model only)
Power	230V AC, 50Hz

Models, Features and Specifications subject to change without notice.