



Bio Optics

P R O D U C T C A T A L O G

· New runners of Medical market ·

BIO Plug™

Dissolvable Lacrimal Punctum Plug




10 pcs/box



BIO Plug

- Use to Polydioxanone Material
- Violet Color
- Never touch the eye
- Comfortable after proper placement
- Absorption of plug approximately 180days
- Available in two size: 0.4mm and 0.5mm diameter



Indication

- Dry eye
- Refractive surgery
- Contact lens wear
- Seasonal allergy
- Cataract Surgery



BIO Intubation Set™

Lacrimal Intubation Set




5 sets/box

BIO Lacrimal Intubation Set


- New state art of Lacrimal Intubation system.
- Bio compatible silicone tubing.

Stainless steel probes with silicone tubing



BS-1001, BS-2001

Probes - 11 cm length X 0.40 mm diameter (4 1/4 in X 27G)
Silicone tubing - 30 cm length, 0.30 mm inside diameter,
0.64mm outside diameter



BS-1002, BS2002

Probes - 11 cm length X 0.60 mm diameter (23G x 7 in)
Silicone tubing - 30 cm length, 0.30 mm inside diameter,
0.64 mm outside diameter



· New runners of Medical market ·

B10 Ophthalmic Micro Surgery Knife

- New Feature of Laser marking Depth Indicator in Keratome Knives.
- Better to Accurate Cut : Laser point Mark

LASER MARKINGS

1.50mm
1.75mm
2.00mm



BIO OPHTHALMIC MICRO KNIFE

01. BIO Keratome Knife



Model No.	Size	Description
B36-2200	2.20mm	Phaco Slit Knife, Angled, Double-Bevel
B36-2600	2.60mm	
B36-2750	2.75mm	
B36-2800	2.80mm	
B36-3000	3.0mm	

02. BIO MVR Knife



Model No.	Size	Description
B36-1900	1.1mm	19G MVR Knife, Straight
B36-2000	0.9mm	20G MVR Knife, Straight

03. BIO Crescent Knife



Model No.	Size	Description
B36-4805	2.0mm	Bevel-up Crescent Knife

04. BIO Scleral Knife



Model No.	Size	Description
B36-5700	2.0mm	200 / 3.2

05. BIO Side Port Knife



Model No.	Size	Description
B36-7513	1.5mm Depth	15° Stab Knife, Straight

BIO Color Tear Test™

Schirmer Strip



100 pcs/box

- Easy to use
- Blue-black dye line
- Different shapes for right(R) and left(L) eye




Critical Value

- Normal which is 10mm-30mm wetting of the paper after 5minute.
- Abnormal which is \leftarrow 10mm wetting of the paper after 5minute.

BIO PVA Eye Spear™

- PVA Eye Spear
- Stiff without water and soft with water
- Provide high safety to the cornea
- Micro point absorbent
- Lint free
- Plastic handle



100 pcs/box

BIO Micro Syringe Filter™

Bio PES Syringe Filter

Specification

- 0.2mm x 25mm
- Safety to PH1~14
- Hydrate Membrane Filter
- Gamma Radiation



50 pcs/box



· *New runners of Medical market* ·

Bio Optics는 최상의 제품생산, 최상의 고객만족
최상의 서비스를 지향합니다.

Bio Optics Co., Ltd.

No.1003, SICOX Tower,484 Dunchon-daero, Jungwon-gu,
Seongnam-si, Gyeonggi-do, KOREA 13229

T 82-31-747-4850 F 82-31-747-4820

www.biooptics.co.kr