

# HERASHAPE



# HERASHAPE

HERASHAPE is a non-invasive laser procedure to destroy fat cells in the abdomen and the flank by using 42 ~ 47°C of Heat.



Non-invasive treatment

Non-surgery treatment

Quick return to normal life

Short treatment time (all 25 minutes)

A little pain.

# HERASHAPE


- Non-invasive treatment.
- Non-surgery treatment.
- No downtime without necessity of recovery time.
- In 25 minutes for one time of treatment
- Feel just a little tingling sensation



Laser Diode : MADE IN GERMANY

# HERASHAPE

The radiated laser energy by **1060nm** and **635nm** raises temperature up to **42 ~ 47°C** and stimulates fat cells. Stimulating disintegrate fat cells and induce apoptosis on the basis of the principle of hyperthermia therapy.



**Target**


The 1060nm,635nm wavelength selectively reaches subcutaneous fat cells below the dermal layer

**Destruction**


Heat(42~47°C) destroys only fat cells

**Removal**


Destroyed fat cells are removed naturally in the body




Choice the area of body for treatment among abdomen, flanks, forearm, thigh etc.



Apply 42~47°C heat to subcutaneous fat layer and induce apoptosis.



Fat cells are gradually removed from body by natural metabolism.



Fat layer thickness decrease continuously over a few weeks.

# HERASHAPE



- With cooling system in HERASHAPE, it prevents burning on skin surface.
- The laser is only emitted when the handpiece is safely contacted to four contact sensors to protect skin from burning.
- The 635nm wavelength of LLLT(Low Level Laser Therapy) works simultaneously to reduce fat together with 1060 laser. LLLT biologically decays fat.
- If device is in abnormal status while operation, press emergency button then all operation will be stopped.
- If a patient feels pain during treatment, press nursing call button in person then laser emitting will be stopped and cooling system will be running, which will make a patient comfortable for one minute.



## BEFORE & AFTER



Minimize procedure time. 25 minutes



several parts can be treated at the same time.



Effect during 8 to 12 weeks.



Non-surgery treatment

# HERASHAPE



**Elimtek Co., Ltd.**

R401, SKN Techno Park Bizcenter, 124,  
Sagimakgol-ro Jungwon-gu, Seongnm-si Gyonggi-do, Republic of Korea

T. +82-31-742-7128 F. +82-31-744-2627 E. [info@elimtek.co.kr](mailto:info@elimtek.co.kr)

[www.elimtek.co.kr](http://www.elimtek.co.kr)