



Supreme

Nurse Call System




Innovative and integrated
Nurse Call Solution for
Hospitals and Healthcare centers



SUPREME Nurse Call System is a nursing management system that can be used for both wired and wireless connections. The features of this system make it ideal for use in

Increases the safety and comfort of patients and residents. The system is comprised of a set of push buttons and wireless sensors which make it possible for nursing staff to keep track at all times of the state of each patient or resident, thereby increasing their



Wireless System

The wireless feature provides flexibility, expandability and ease of installation. and most importantly it is the future of the NCS enables the system to be used in



Wired /Wireless IP based nurse call system



Supreme

Nurse Call System..

Easy To Use

Quick & Easy Installation


Safe & Secure

Cost Effective

Nurse Station/Display Unit : The Nurse Call System or NCS eradicates the frustration of patients to call the nurses. NCS enables even small hospitals and nursing homes to have dedicated system to call nurses. NCS is economical for hospitals and nursing homes. The system eases the management of patient related services like Blue Code and



Ds1



Ds2



Ds3




Ds4

Bed Unit & Handset : The Bed Unit is located near the patient bed. The unit has CALL, WAIT, BLUE CODE and RESET buttons and female socket to attach the Handset. The Bed Side Unit has battery for wireless models.



Foot operated
Fb1




Wall mounted Call switch
Cb1



Wearable pendant
Wb1



Bed side unit
Bm1



Bed side unit
Bs1

Information is sent to fixed control posts, computers and wireless terminals. This enables communication between patients and nurses regardless of location, thus optimizing the work of the nursing staff.



Call Logger

Date	Time	Ward number	Bed number	Status
2/2/2016	22:11:25	A-1	101	call
2/2/2016	22:12:17	A-3	117	call
2/2/2016	22:13:20	A-1	101	Cancel
2/2/2016	22:15:17	A-2	108	call
2/2/2016	22:12:19	A-3	117	Cancel

The Call Logger will automatically log and store all nurse call logs on board. The data can be viewed on any computer. The date, time and room location are stored and easily searched. The data file can be viewed by removing the sd card and inserting it in any device, PC, Mac, Tablet or

Door Unit : The door unit indicates the patient call in the corridor. The red color of unit indicates patient call.

Toilet Pull Cord : The NCS has toilet pull cord in toilets of hospital for emergency situations like heart attack of

DEALER STAMP

WWW.WIRELESS4U.IN
WWW.SIWORLD4U.COM

SUPREME INTERNATIONAL
4772-4773, T-1, 3RD FLOOR,
BHARAT RAM ROAD,
23 DARYA GANJ, NEW DELHI