

WATER BATH PARAFFIN (LI-WBPF-14B)

APPLICATION :

The Paraffin Section Mounting Bath handles paraffin wax samples in pathology and histology labs. The hot distilled water floating out bath lets user meticulously manipulate and locate paraffin wax section on glass slide(s). Excessive wax or damage to sections can be removed with a filter paper on the surface of the water.



FEATURES :

- Bath made of SS. Outer made of MS duly powder coated.
- Thermostat with/or Energy regulator maintains the temperature of the bath.
- Double walled with cups and sleeves for glass tubes.
- Industrial grade energy efficient heating elements reduce power bill drastically.
- Drain nozzle provided.

TECHNICAL SPECIFICATIONS :

Model : LI-WBPF-14B			
Cups & Sleeves	4 cups & 6 Sleeves	6 cups & 8 sleeves	12 cups & 14 sleeves
Constructions	Bath of SS and outer body of MS duly powder coated		
Temp. Controller	Thermostatic or Digital Controller (optional)		
Temp Range	Above ambient 5°C to 95°C		
Temp. Accuracy	±1°C		
Power	230V AC, 50Hz		

OPTIONAL :

- Digital Temperature Indicator-cum-Controller.

Models, Features and Specifications subject to change without notice.