



# ANILINE POINT APPARATUS-METHOD A-ASTM D611, IP 2, IS 1448-HAND OPERATED

## Compliance with following International Standards

ASTM D611, IP 2

### Scope

These test methods cover the determination of the aniline point of petroleum products and hydrocarbon solvents.

Method A is suitable for transparent samples with an initial boiling point above room temperature and where the aniline point is below the bubble point and above the solidification point of the aniline sample mixture.

### Brief Construction Details

Consists of following Parts and accessories.

#### 1) Test Tube

Diameter : 25 mm  
Length : 150 mm  
Material : Made of heat-resistant glass.



#### 2) Jacket

Diameter : 37 - 42 mm in diameter  
Length : 175 mm  
Material : Made of heat-resistant glass

#### 3) Stirrer

Manually operated  
Made of Stainless steel material  
Diameter : 2 mm

#### 4) Cork for tube

#### 5) Stand support

### Optional Accessories at extra cost:

Spare outer Jacket with cork @ Rs 825.00  
Spare tube with cork @ Rs 255.00

#### Corporate Office:

A-1301, BVR Ek, Opp. Hotel Inder Residency,  
Near Gujarat College, Ellisbridge,  
Ahmedabad-380 006.  
Phone : (079) 6621 1234  
Website : [www.eieinstruments.com](http://www.eieinstruments.com)

#### Calibration Laboratory:

B-14, Zaveri Industrial Estate, Opp. Shyam Villa,  
Singarva-Kathwada Road,  
Kathwada, Ahmedabad-382 430.  
Phone : (079) 6604 0660 • Fax : (079) 6604 0600  
E-mail : [info@eieinstruments.com](mailto:info@eieinstruments.com)

GST TIN No. : 24AABCE4018L1ZU  
I.T. PAN No. : AABCE 4018L  
IEC No. : 0813026598 Dt. 04/03/2014  
CIN No. : U29199GJ2004PTC045078  
TAN No. : AHME00572F



# EIE INSTRUMENTS PVT. LTD.

(An ISO 9001:2015 Certified Company)

IP 20C (-38/42°C x 0.2C)	@ Rs 850.00
IP 21C (25/105°C x 0.2C)	@ Rs 850.00
IP 59C (90/170°C x 0.2C)	@ Rs 850.00

#### Corporate Office:

A-1301, BVR Ek, Opp. Hotel Inder Residency,  
Near Gujarat College, Ellisbridge,  
Ahmedabad-380 006.  
Phone : (079) 6621 1234  
Website : [www.eieinstruments.com](http://www.eieinstruments.com)

#### Calibration Laboratory:

B-14, Zaveri Industrial Estate, Opp. Shyam Villa,  
Singarva-Kathwada Road,  
Kathwada, Ahmedabad-382 430.  
Phone : (079) 6604 0660 • Fax : (079) 6604 0600  
E-mail : [info@eieinstruments.com](mailto:info@eieinstruments.com)

GST TIN No. : 24AABCE4018L1ZU  
I.T. PAN No. : AABCE 4018L  
IEC No. : 0813026598 Dt. 04/03/2014  
CIN No. : U29199GJ2004PTC045078  
TAN No. : AHME00572F