



# Gel Card Centrifuge & Incubator



**OmniCard 12 C**  
*Efficiently Flexible*



**OmniCard 12 i**  
*Efficiently Flexible*



**Universal Gel Card System**  
*Provides flexibility to run all available gel cards*



**REMI Omnicard** is an efficiently flexible gel card system provides you freedom to select any brand of gel card with replicating standard results. It provides you the various option of run parameters for centrifugation & incubation. **Omnicard** is an open gel card system for economically viable solution for blood grouping, Coombs tests (DAT) & cross matching, etc.

**Omnicard 12 C** centrifuge fits all types of gel cards, with flexible 10 program setting (speed & time) including one sequential program.

**Omnicard 12 i** incubator fits all types of gel cards with flexible temp. & time setting. Additional dry bath slots for 5-7 ml test tubes makes it helpful for tube agglutination methods also.



Specifications subject to change without notice



# Gel Card Centrifuge & Incubator

## Salient features

### OmniCard 12 C

- Centrifugation of 12 Gel cards.
- Digital user interface screen provided for control and visualization of rpm and time.
- Maximum speed 1600 rpm.
- Maximum RCF 279g±1%
- 8 open program for setting RPM, acceleration & deceleration, time with additional two sequential sets of parameters
- Audio-Visual alarm provided for Imbalance & end of cycle.
- Universal for all types of standard gel cards.
- Speed control through VFD and dynamic brake system.
- Light weight yet robust design.
- Ergonomically designed shape for ease of operation
- Single Hand Rotor removal for cleaning

### OmniCard 12 i

- Incubation of 12 Gel cards & 12 test tubes simultaneously.
- Digital user interface screen provided for control and visualization of remaining incubation time and real temperature.
- Fixed Pre-set temperature 37°C (Pre heating within 8 minutes)
- Fixed pre-set incubation time of 15 minutes and can be adjustable from 0 to 60 min.
- Audio-Visual alarm provided for High/Low temperature & End of Cycle.
- Auto Pop-up Card mechanism, when lid opened.
- Light weight yet robust design.

## Technical Specifications

Parameter	OmniCard 12 C Gel Card Centrifuge	OmniCard 12 i Gel Card Centrifuge
Programs	Adjustable RPM & Time as specified by Gel card Supplier	Adjustable incubation Temperature & Time as specified by Gel card Supplier
Run specification	Max Speed : 1600±10 RPM Max RCF : 279±1% g	37°C ±1°C Adjustable time up to 60 minutes
Net Weight	18 Kg.	4 Kg.
L x D x H	420 x 320 x 286 mm	300 x 250 x 150 mm
Audio Visual Alarms	Imbalance & End of cycle	High-low temperature, End of cycle
Capacity	12 Gel cards	
User interface	LCD Display with Touch key pad	
Electrical	220-240 volts, 50 Hz, SINGLE PHASE	



Specifications subject to change without notice