



Dual °C & °F



Fever Alarm



Accurate Result



Mercury Free

DT-01KL Digital Thermometer

Features:

- Non Mercury (100% Safe).
- Clinically Tested.
- Highly Accurate.
- Alarm Signal and Memory Function.
- Measurement in °C & °F.
- Easy to Read, LCD Display.
- Can be Used in Oral Cavity, Rectal & Under the Armpit.

Specifications:

Display Type	LCD
Battery	DC 1.5V Type-LR41
Resolution	0.1
Storage Conditions	-25°C -55°C ≤85% RH
Measuring Range	35.5°C - 42°C /(95.9-109.4)°F
Max Error	0.1 °C/0.2°F
Working Conditions	10°C - 35°C ≤85%RH

How to Use:



Press Power Button to turn on the unit.



Place the thermometer under the tongue or the underarm & And wait for beeper sound.



Read Temperature.



Press Power Button to turn off the unit or leave it aside it will automatically off after some time.

Contact us:



Mobile No./ WhatsApp No. : +91 9310113243



Web : www.klifecare.com



Email : info@klifecare.com

Kannu Impex (India) Pvt. Ltd.

412, Suncity Trade Tower, Sector-21, Gurugram-122015