

DIGITAL COLONY COUNTER (LIM-163 & LIM-165)

APPLICATIONS :

Ideal for all microbiology applications, it provides fast and accurate counting of bacterial and mould colonies.

This colony counter is easy to use. Simply place the petri-dish on the glow plate and touch the dish with a felt tip pen to mark each colony. The touch pressure causes a count to be registered on the digital display and an audible tone confirms each count made.



TECHNICAL SPECIFICATIONS :

| | Model : LIM-163 | Model : LIM-165 |
|----------------|---------------------------------------|--------------------------------------|
| Display | 4 digit LED display (0 to 9999 count) | 3 digit LED display (0 to 999 count) |
| Glow dish Size | 110mm | |
| Light source | 25W (Milky Candle Lamp) | |
| Power | 230V (±10%) AC, 50Hz | |

Models, Features and Specifications subject to change without notice.