

Company Information

Company Name SAMMYUNG TECH CO.,LTD

| | | | |
|---------|----------------------|----------------|------------------|
| CEO | KANG TAE-SEOK | Position | CEO |
| | | Contact Person | KANG TAE-SEOK |
| Website | www.sammyungtech.com | Tel | +82-31-498-8550 |
| Tel | +82-31-498-8550 | E-mail | ank7772naver.com |
| E-mail | ank777@naver.com | | |

● Company Introduction

Sammyung Tech Ltd., is a beauty medical device manufacturer established in August 2014 and is a venture company with a technical development department and have 3 patent registrations and 2 patent applications. The company obtained ISO9001, 14001, ISO13485 certification. It was selected as an export frontier company in 2015 and 2017 and has been designated as a Gyeonggi export business in 2017 and has been designated as a Gyeonggi business in 2018. Our main products include low-frequency foot massage devices, wearable smart belts, and skin beauty device; which provide optimal skin care as they have temperature measurement function and hot massage cooling massage function.

● Items (품목 및 제품명 등 한줄로 간단한 소개)

- 1)Low frequency massager 2)Low frequencysmart belt
- 3)Multi-skinned beauty device

● Features

(Item 1)
Low frequency massager



(photo)

(Item 2)
Low frequency smart belt



(photo)

● Specification

This product is effective in relieving pain by spreading fine low frequency currents through the human body to promote muscle movement and blood circulation. Its main features include the ability to use foot massagers and body pads at the same time, eight types of auto mode, 18 types of manual mode, and 99 levels of fine-tuning.

The wearable low-frequency smart belt measures the amount of exercise and calorie consumption depending on the time and intensity of the device. And you can compare, record, and save the result data using the app on your smartphone.

(Item 3)
Multifunctional skin
beauty
device



- 1) cool massage helps to close widened pores sweating out in summer.
- 2) hot massage helps blood circulation.
- 3) led lights anti-wrinkle using LED lights.
- 4) Micro-current have young and healthy skin with collagen.
- 5) vibration don't tap with your fingertips.
- 6) Temperature measurement
What would the temperature of my skin?